

# Generation of DNA fragments by DNase digestion

## William H. Wilson and Declan Schroeder

# **Abstract**

Generation of DNA fragments by DNase digestion involves digesting DNA for a range of times, then picking the time that gives optimal-sized DNA fragments (typically 1000–4000 bp).

Citation: William H. Wilson and Declan Schroeder Generation of DNA fragments by DNase digestion. protocols.io

dx.doi.org/10.17504/protocols.io.dy47yv

Published: 18 Jan 2016

## **Protocol**

## Step 1.

In a 50  $\mu$ L reaction volume, resuspend 8  $\mu$ g DNA in 50 mM Tris·HCl (pH 7.6), 10 mM MnCl<sub>2</sub>, 100  $\mu$ g mL<sup>-1</sup> bovine serum albumin, and 0.01 SU mL<sup>-1</sup> DNase I.

## Step 2.

Remove 5  $\mu$ L aliquots (adding to 45  $\mu$ L TE buffer, pH 7.6) 0, 0.5, 1, 2, 5, 10, 15, and 30 min after addition of the digestion mixture.

## Step 3.

Immediately transfer to a tube containing 25 µL Tris-buffered (pH 7.0) phenol.

### NOTES

#### Declan Schroeder 12 Oct 2015

Typically the shorter incubations, up to 2 min, give optimally sized fragments.

#### Step 4.

Perform a phenol:chloroform (1:1) extraction.

#### Step 5.

Perform a chloroform extraction.

#### Step 6.

Perform a chloroform extraction once more.

#### Step 7.

Precipitate the fragmented DNA.

#### Step 8.

Wash with 70% ethanol.

# Step 9.

Dry.

#### **Step 10.**

Resuspend fragmented DNA in 23  $\mu$ L of Blunt-ending Mix (100  $\mu$ M dNTPs, 1  $\times$  T4 DNA Pol Buffer).

#### Step 11.

Heat at 65°C for 30 min to resuspend DNA and inactivate any DNase I that was carried over.

© DURATION

00:30:00

## Step 12.

Cool to room temperature.

## **Step 13.**

Add 2.5 U Klenow fragment and 5 U T4 DNA polymerase.

## **Step 14.**

Incubate the reaction at 37°C for 1 h.

**O DURATION** 

01:00:00

## **Step 15.**

The fragmented and blunt-ended virus DNA can be run on a 1% agarose gel prior to excising fragments in the 1000–4000 bp range using a standard gel extraction procedures before downstream cloning.

## **P** NOTES

# **Declan Schroeder** 12 Oct 2015

NB do not excise fragments smaller than 1000 kb, as downstream cloning will preferentially clone the smaller fragments.