# **Fungal DNA Lysis Buffer**

#### **Adam Taranto**

### **Abstract**

Lysis buffer for extraction of DNA from fungal material. Contains ionic and non-ionic detergents. Contains RNAse-A.

### Composition:

NaCl 150mM EDTA pH 8.0 1 mM Tris Base pH 8.0 10 mM Polyvinylpolypyrrolidone 2% (PVPP) IGEPAL-CA630 0.5% Sodium deoxycholate 0.5% **SDS** 1% RNAse-A 100 μg/mL

Citation: Adam Taranto Fungal DNA Lysis Buffer. protocols.io

dx.doi.org/10.17504/protocols.io.fvfbn3n

Published: 17 Oct 2016

# **Materials**

0.5mm diameter glass beads <u>SI-BG05</u> by <u>Scientific Industries, Inc.</u> IGEPAL-CA630 <u>I3021 SIGMA-ALDRICH</u> by <u>Sigma Aldrich</u>

- ✓ Tris(Hydroxymethyl)aminomethane PubChem CID: 6503 by Contributed by users
- ✓ Sodium Dodecyl Sulfate <u>PubChem CID: 3423265</u> by Contributed by users
- ✓ EDTA Disodium Salt <u>PubChem CID</u>: 8759 by Contributed by users
- ✓ Sodium Chloride <u>PubChem CID: 5234</u> by Contributed by users
- ✓ Sodium Deoxycholate <u>PubChem CID: 23668196</u> by Contributed by users RNase A (10 mg/mL) <u>EN0531</u> by <u>Thermo Fisher Scientific</u>
- Polyvinylpolypyrrolidone CID 6917 by Contributed by users

# **Protocol**

#### **Prepare Stock Solutions**

Step 1.

10% SDS solution

Dissolve 10 g SDS in 80 mL ddH<sub>2</sub>O Adjust volume to 100 mL

**■** AMOUNT

10 g Additional info:



**REAGENTS** 

✓ Sodium Dodecyl Sulfate <u>PubChem CID: 3423265</u> by Contributed by users

### **Prepare Stock Solutions**

Step 2.

### 1M Tris Base, pH 8.0

Dissolve 12.114g Tris-Base in 80 mL ddH<sub>2</sub>O

pH to 8.0 with 10 M HCl Adjust volume to 100mL



**■** AMOUNT

12 g Additional info:



✓ Tris Hydroxymethylaminomethane PubChem CID: 6503 by Contributed by users

### **Prepare Stock Solutions**

# Step 3.

### 0.5 M EDTA, pH 8.0

Add 18.61 g of disodium EDTA • 2H<sub>2</sub>O to 80 mL of ddH<sub>2</sub>O Adjust the pH to 8.0 with NaOH (20 NaOH pellets) Adjust volume to 100 mL



19 g Additional info:



**REAGENTS** 

✓ EDTA Disodium Salt <u>PubChem CID</u>: 8759 by Contributed by users.

### **Prepare Stock Solutions**

# Step 4.

### 3 M NaCl

Dissolve 17.532 g NaCl in 80 mL ddH<sub>2</sub>O

Adjust volume to 100 mL



18 g Additional info:



✓ Sodium Chloride <u>PubChem CID</u>: 5234 by Contributed by users

# Prepare Lysis Solution - 100 mL

# Step 5.

Add 500µL IGEPAL-CA630 to 82.3 mL ddH<sub>2</sub>O

Stir until dissolved

AMOUNT

500 µl Additional info:



**REAGENTS** 

IGEPAL-CA630 <u>I3021 SIGMA-ALDRICH</u> by <u>Sigma Aldrich</u>

NOTES

Adam Taranto 22 Sep 2016

Final concentration: 0.5%

### Prepare Lysis Solution - 100 mL

Step 6.

Add 0.5g Sodium deoxycholate

Stir until dissolved



REAGENTS

✓ Sodium Deoxycholate PubChem CID: 23668196 by Contributed by users

NOTES

Adam Taranto 22 Sep 2016

Final concentration: 0.5%

# Prepare Lysis Solution - 100 mL

Step 7.

#### **Add PVPP**

Add 2g Polyvinylpolypyrrolidone, stir until dissolved.



2 g Additional info:



✓ Polyvinylpolypyrrolidone <u>CID 6917</u> by Contributed by users

NOTES

Adam Taranto 17 Oct 2016

Final concentration 2%

Adam Taranto 17 Oct 2016

PVPP precipitates phenols, tannins, alkaloids which may inhibit enzymes or damage DNA

# Prepare Lysis Solution - 100 mL

Step 8.

Add 10 mL of 10% SDS stock solution

### Stir until dissolved

**■** AMOUNT

10 ml Additional info:

REAGENTS

✓ Sodium Dodecyl Sulfate PubChem CID: 3423265 by Contributed by users

NOTES

Adam Taranto 22 Sep 2016

Final concentration: 1%

Adam Taranto 17 Oct 2016

1% SDS = 34.67 mM

# Prepare Lysis Solution - 100 mL

Step 9.

Add 5 mL of 3M NaCl stock solution

**■** AMOUNT

5 ml Additional info:

REAGENTS

✓ Sodium Chloride PubChem CID: 5234 by Contributed by users

NOTES

Adam Taranto 22 Sep 2016

Final concentration: 150 mM

# Prepare Lysis Solution - 100 mL

Step 10.

Add 200 µL of 0.5 M EDTA (pH 8.0) stock solution

**■** AMOUNT

200 µl Additional info:

REAGENTS

✓ EDTA Disodium Salt <u>PubChem CID</u>: 8759 by Contributed by users

NOTES

Adam Taranto 22 Sep 2016

Final concentration: 1 mM

# Prepare Lysis Solution - 100 mL

### **Step 11.**

Add 1 mL Tris-Base (pH 8.0)

Stir until dissolved

**■** AMOUNT

1 ml Additional info:



✓ Tris Hydroxymethylaminomethane PubChem CID: 6503 by Contributed by users

NOTES

Adam Taranto 22 Sep 2016

Final concentration: 10 mM

# Prepare Lysis Solution - 100 mL

**Step 12.** 

Add 1mL of RNase A (10 mg/mL) stock solution

**■** AMOUNT

1 ml Additional info:



RNase A (10 mg/mL) EN0531 by Thermo Fisher Scientific

NOTES

Adam Taranto 22 Sep 2016

Final concentration: 100 µg/mL

# Prepare Lysis Solution - 100 mL

**Step 13.** 

Store lysis solution at 4°C