



Jul 24, 2019

## Prepare NGM plates for nematode, with peptone, without fungizone V.2

Cancer Research UK / Wellcome Gurdon Institute media kitchen<sup>1</sup><sup>1</sup>Wellcome Trust / Cancer Research UK Gurdon Institute

1

Works for me

[dx.doi.org/10.17504/protocols.io.5smg6c6](https://doi.org/10.17504/protocols.io.5smg6c6)

Cristian Riccio ⚡

### ABSTRACT

Prepare NGM plates for nematodes, with peptone, without fungizone

### MATERIALS

NAME ▾

CATALOG # ▾

VENDOR ▾

double distilled water (ddH<sub>2</sub>O)

Sodium chloride meets analytical specification of Ph.Eur Fisher Chemical

S/3160/65

Fisher Scientific

SYCHEM autoclave

View

Syschem

Agar

AGA02

Formedium

Minisart filters pore size 0.2 µm

16534K

Sigma Aldrich

Bacto™ Peptone

211677

Thermo Fisher Scientific

### SAFETY WARNINGS

Make sure you know how to use the autoclave before starting this protocol.

### 1 NGM\_agar\_1L.xls

### 2

Ingredients		Quantity
NaCl		3g
bacto peptone		2.5g
Agar		17g
Double distilled H <sub>2</sub> O	972ml	

### 3

Measure 972ml double distilled H <sub>2</sub> O and put in a 1L duran bottle with a magnetic flea
Add 3g NaCl, 2.5g bacto peptone and stir.
Stir until all solutes are dissolved then add the agar.
Leave magnetic flea in bottle.
Label, date and autoclave



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited