

# Nuclei Isolation from Human Brain Using Sucrose Gradient

Fatma Ayhan, Genevieve Konopka

## Abstract

This protocol outlines our preparation of single-nuclei suspension from surgically acquired fresh human adult brain tissue.

**Citation:** Fatma Ayhan, Genevieve Konopka Nuclei Isolation from Human Brain Using Sucrose Gradient. **protocols.io**  
dx.doi.org/10.17504/protocols.io.scneave

**Published:** 03 Aug 2018

## Protocol

### Step 1.

Prepare Nuclei Suspension Buffer (NSB)

For drop-seq: 1x PBS, 0.01% BSA (ultra pure), 0.2U/ul RNase inhibitor.

For 10X Genomics: 1X PBS, 1%BSA (ultra pure), 0.2U/ul RNase inhibitor

### Step 2.

Prepare lysis buffer and sucrose solution:

Lysis buffer	30 ml	Final
2M Sucrose	4.8 ml	0.32 M
1M CaCl <sub>2</sub>	150 µl	5 mM
1M Mg(Ac) <sub>2</sub>	90 µl	3 mM
0.5M EDTA	6 µl	0.1 mM
1M Tris-HCl(pH8.0)	300 µl	10 mM
0.1M PMSF	30 µL	0.1 mM
100% Triton X-100	30 µL	0.10%
100% NP-40	30 µL	0.10%
Protease inhibitor	300 µL	1X
RNase inhibitor	150 µL	0.1 U/µL
Molecular biology grade (MBG) water	24.26 ml	

<b>Sucrose solution</b>	<b>50 ml</b>	<b>Final</b>
2M sucrose	45 ml	1.8 M
1M Mg(Ac) <sub>2</sub>	150 µl	3 mM
1M Tris-HCl(pH8.0)	500 µl	10 mM
MBG water	4.35 ml	

### **Step 3.**

Put 250mg of brain tissue into a douncer with 10ml of ice-cold lysis buffer on ice.

### **Step 4.**

Dounce the tissue for 2 min while on ice.

### **Step 5.**

Take 18ml of sucrose solution into a new 30ml clear ultracentrifuge tube on ice.

### **Step 6.**

Place the homogenized brain into ultracentrifuge tube on ice (gradient with the homogenized brain on top of the sucrose solution).

### **Step 7.**

Weigh the ultracentrifuge tubes, and adjust the weight by lysis buffer.

### **Step 8.**

Ultracentrifuge with SW28 rotor at 24,400RPM (=107,163.6 RCF) for 2.5h at 4°C.

### **Step 9.**

Remove supernatant, add 500µL of NSB to each pellet, and incubate them on ice for 10min.

### **Step 10.**

Resuspend the nuclei in NSB, filter through Flowmi cell strainer (Bel-Art, H13680-0040) and transfer the suspension into one 2ml tube on ice.

**Step 11.**

Stain 10 $\mu$ l of nuclei suspension with DAPI (1:5000) or 0.4% Trypan Blue to count.

**Step 12.**

Adjust the volume with NSB to 1000 nuclei/ $\mu$ l (for 10X Genomics) and 300 nuclei/ $\mu$ l (for drop-seq).