



Feb 01, 2020

RNA Extraction for RIN and DV 200 Analysis

Jamie Allen¹, Elizabeth Neumann¹, Maya Brewer¹, Jeff Spraggins¹, Danielle Gutierrez¹, Mark de Caestecker²

¹Vanderbilt University, ²Division of Nephrology, Vanderbilt University Medical Center

1 Works for me dx.doi.org/10.17504/protocols.io.86nhzde

VU Biomolecular Multimodal Imaging Center

Human BioMolecular Atlas Program (HuBMAP) Method Development Community



ABSTRACT

Scope: Extract RNA for RIN and DV 200 Assessment.

Expected Outcome: RIN and DV 200 measurements for tissue quality assessment.

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Adapted from: Masato Hoshi, MD, PhD, Jain Lab, Division of Nephrology, Washington University School of Medicine

MATERIALS TEXT

Materials

1. Forceps
2. 1.5mL microcentrifuge tubes
3. 23G and 27G needle syringes
4. TRIzolLS
5. DEPC treated water
6. RNase free water
7. Phase Lock Gel, 5PRIME #2302830
8. Chloroform
9. Centrifuge
10. Isopropanol
11. Ethanol
12. Qiagen RNeasy Plus Micro Kit
13. 2-ME
14. RPE Buffer
15. Bioanalyzer

Solutions:

1. Glycogen Solution
5 mg/mL

Protocols:

Qiagen RNeasy Plus Micro Kit

Homogenize Tissue

- 1 Place tissue section in 750 μ L of TRIzolLS.

- 2 Homogenize the sample:
Use a homogenizer or pass sample through a 23G then a 27G needle syringe.
- 3 Freeze samples on dry ice and re-thaw once. Mix well by vortexing.
- 4 Add RNase free water to a final volume of 960 μ L.
- 5 Add 40 μ L of glycogen solution (5mg/mL, making a final conc of 200 μ g/mL of glycogen) to each tube and mix well.


Loading

- 6 Incubate for 5 minutes at 20°C (room temperature).
- 7 Add 250 μ L of chloroform and shake for 15 seconds.
- 8 Centrifuge at 12,000 x g for 10 minutes at 4°C.
- 9 Move supernatant to a new 1.5 mL tube.
- 10 Add 600 μ L isopropanol.
- 11 Vortex and incubate the tube at -20°C for 20 minutes.
- 12 Centrifuge at 20,000 x g for 20 minutes at 4°C.
- 13 Discard supernatant and add 600 μ L 80% Ethanol.
- 14 Vortex and centrifuge at 20,000 x g for 5 minutes at 4°C
- 15 Discard supernatant and air dry for a few minutes.
- 16 Add 300 μ L of RLT plus buffer (from Qiagen RNeasy Plus Micro Kit) containing 2-ME (10 μ L of 2-ME for 1 mL of RLT Plus Buffer, 1% 2-ME) and mix.
- 17 Add 450 μ L 100% ethanol (1.5 volumes) and mix well.
- 18 Transfer 750 μ L of the sample to RNeasy MinElute spin column and centrifuge for 1 min at 12,000 x g at 24°C.

Collection

- 19 Centrifuge Phase Lock Gel to move gel to the bottom of column (12,000 x g for 5 minutes).
Then transfer the TRIzol mixture mixture to the column.
- 20 Add 500 μ L of RPE buffer to the column and centrifuge at 12,000 x g for 1 min at 24°C.
- 21 Discard flow through and add 500 μ L RPE buffer to the column. Centrifuge at 12,000 x g for 2 min at 24°C.
- 22 Transfer column to a new collection tube. Discard the old collection tube with the flow through.
- 23 Transfer the collected flow through back through the same column and centrifuge for 1 min at 12,000 x g at 24°C. (This is a double application of sample). This time discard the flow.

- 24 Open the lid of the column and centrifuge at 12,000 x g for 5 min at 24°C to dry up the column membrane completely.
- 25 Place the column in a new 1.5 mL tube and add 10 uL of RNase free water on the center of the column membrane. Incubate for 2 min at 20°C.
- 26 Centrifuge at 12,000 x g for 3 min at 24°C.
- 27 Add 10 uL of RNase free water again on the center of the column membrane and incubate for 2 min at 20°C.
- 28 Centrifuge at 12,000 x g for 3 min at 24°C (this is a double elution).
- 29 Total volume should be ~20 uL.
- 30 Provide VANTAGE Core with 10 uL of solution for RIN and DV 200 analysis.

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited