



Sep 17, 2019

## Genomic DNA extraction from diatom *P. multistriata* V.1

FRANCESCO MANFELLOTTO<sup>1</sup>, monia teresa russo<sup>1</sup>, Rossella Annunziata<sup>2</sup>, Mariella Ferrante<sup>2</sup>

<sup>1</sup>SZN, <sup>2</sup>Stazione Zoologica Anton Dohrn

1 Works for me dx.doi.org/10.17504/protocols.io.7dghi3w



FRANCESCO MANFELLOTTO  
SZN




### ABSTRACT

Genomic DNA extraction from diatom *P. multistriata*

- 1 Grow cells in 250 ml
- 2 filter cultures with 1.2 µm nitrocellulose membranes
- 3 Centrifuge at 6000 rpm for 5 minutes at 4 °C and remove medium
- 4 Resuspend cells with 500 µl of TE buffer (10 mM TrisHCl pH 7.6 and 1 mM EDTA pH 8.0)
- 5 Add:
  - 400 mg of 0.2-0.3 mm zirconia/silica diameter beads,
  - 500 µl phenol.
- 6 rinse filter with 1 ml f/2 medium in a falcon and then move cells to a 2 ml eppendorf
- 7 Mix with vortex 30 Hz 3 times for 85 seconds, each time put sample in ice for 60 seconds before vortex.
- 8 · Centrifuge at 10000 rpm for 5 minutes at 4°C.
- 9 Recover aqueous phase in new 1.5 ml eppendorf (about 600 µl)
- 10 Add 500 µl of PCI (Phenol:Chloroform: isoamyl alcohol 25:24:1 v/v) and mix by inversion.

- 11 Centrifuge at 10000 rpm for 5 minutes at 4°C
- 12 Move the aqueous phase in a new eppendorf and add 5 µl of RNase-A 10 mg/ml
- 13 Incubate at 37 °C for 30 minutes. 3m
- 14 Add 500 µl di PCI (Phenol:Chloroform: isoamyl alcohol 25:24:1 v/v) mix by inversion.
- 15 Centrifuge at 10000 rpm for 5 minutes at 4°C.
- 16 Move the aqueous phase in a new 2 ml eppendorf and add:
- 50 µl of 3 M NaAc (pH ± 5)
  - 1 ml of ethanol 96% (- 20 °C)
  - 2 µl glycogen (- 20 °C)
- 17 incubate over night at -20°C. 1h
- 18 Centrifuge the overnight samples at 13000 rpm for 30 minutes at 4°C 3m
- 19 discard aqueous phase
- 20 Wash the pellet by adding 1 ml ethanol 70% and mix gently by inversion
- 21 Centrifuge at 13000 rpm for 10 minutes at 4°C 1m
- 22 discard aqueous phase
- 23 Wash the pellet by adding 1 ml ethanol 70% and mix gently by inversion
- 24 Centrifuge at 13000 rpm for 10 minutes at 4°C 1m
- 25 discard aqueous phase

- 26 Remove aqueous phase and dry pellet at RT for at least 20 minutes 2m
- 27 add 50 µl of Preheated TE 1X (pH 8) or sterile MilliQ water to pellet of DNA  
Preheat the TE pH 8 or sterile MilliQ water at 55 °C
- 28 incubate at 55 °C for 20 minutes 2m
- 29 quantify DNA concentration by nanodrop or Qubit
- 30 in order to check DNA integrity, run a small amount of DNA with 1% agarose gel
- 31 The DNA is ready and store it at -20°C

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited