

# MojoSort™ Positive Selection Columns Protocol

Kelsey Miller

## Abstract

BioLegend MojoSort™ nanobeads work in commonly used separation columns, based on our internal research as well as validation by external testing by academic labs. This simple protocol consists of following the MojoSort™ protocol to label the cells with pre-diluted MojoSort™ reagents and using the columns as indicated by the manufacturer.

Note: Due to the properties of our beads, it may be possible to use far fewer beads than with other commercial suppliers. We recommend a titration to find the best dilution factor. However, as a general rule, dilutions ranging from 1:3 to 1:20 for the Nanobeads can be used.

**Citation:** Kelsey Miller MojoSort™ Positive Selection Columns Protocol. **protocols.io**

dx.doi.org/10.17504/protocols.io.e3kbqkw

**Published:** 06 Jun 2016

## Guidelines

### Commercially available columns:

|                 | Max. number of labeled cells | Max. number of total cells | Cell suspension volume             | Column rinse volume | Cell wash volume | Elution volume |
|-----------------|------------------------------|----------------------------|------------------------------------|---------------------|------------------|----------------|
| Small Capacity  | $10^7$                       | $2 \times 10^8$            | 500 $\mu$ l for up to $10^8$ cells | 500 $\mu$ l         | 3x1 ml           | 1 ml           |
| Medium Capacity | $10^8$                       | $2 \times 10^9$            | 500 $\mu$ l for up to $10^8$ cells | 3 ml                | 3x3 ml           | 5 ml           |
| Large Capacity  | $10^9$                       | $2 \times 10^{10}$         | 5 mL for up to $10^9$ cells        | 20 – 50 ml          | 4x30 ml          | 20 ml          |

### Example of magnetic separation with medium capacity columns:

- Place the column in a magnetic separator that fits the column.
- Wash the column with 3 mL of buffer.
- Add the labeled cell suspension to the column through a 30  $\mu$ m filter and collect the fraction containing the unlabeled cells.
- Wash the column 3 times with 3 mL of buffer and collect the fraction containing the unlabeled cells. Combine with the collected fraction from step c. These cells may be useful as controls, to

monitor purity/yield, or other purposes.

e) Take away the column from the magnet and place it on a tube. Then add 5 mL of buffer and flush out the magnetically labeled fraction with a plunger or supplied device. These are the positively isolated cells of interest; do not discard. To increase the purity of the magnetically labeled fraction repeat the isolation process with a new, freshly prepared column.

**View protocol on Biolegend.com for [Data](#)**

## Protocol

### Step 1.

Prepare a single cell suspension and resuspend the cells with ice cold cell separation buffer (MojoSort™ buffer is recommended).

### Step 2.

Pass the cells through a 70 µm filter, centrifuge (300 x g for 5 minutes), discard the supernatant and resuspend the cells in cell separation buffer. Adjust the cell concentration to  $1 \times 10^8$  cells/mL.

 **DURATION**

00:05:00

### Step 3.

Aliquot 100 µL ( $10^7$  cells) into a new tube.

### Step 4.

Vortex the antibody-conjugated Nanobeads (to resuspend) at max speed, 5 touches, and prepare the dilutions to test. Add 10 µL of pre-diluted conjugated Nanobeads. Mix well and incubate on ice for 15 minutes. Scale up the volume accordingly if separating more cells. For example, add 100 µL of pre-diluted Nanobeads for  $1 \times 10^8$  cells. When working with less than  $10^7$  cells, use indicated volumes for  $10^7$  cells.

Note: Depending on the conjugated nanobead you are using, a wash step may be required here.

 **DURATION**

00:15:00

### Step 5.

Resuspend the cells in appropriate amount of buffer. At least 500 µL is needed for column separation.

Note: There are several types of commercially available columns, depending on your application, choose the one that fits best your experiment. See 'Guidelines' for choosing a column for your

experiment.