

Mar 06,  
2019

Working

## Lipid Biomarker Extraction and Elution into different Fractions from sediment

Lucia Leierer<sup>1</sup>, Antonio V. Herrera-Herrera<sup>1</sup>, Margarita Jambrina-Enríquez<sup>1</sup>

<sup>1</sup>Universidad de La Laguna

[dx.doi.org/10.17504/protocols.io.yg2fye](https://doi.org/10.17504/protocols.io.yg2fye)

Paleochar ERC

 Lucia Leierer  
Universidad de La Laguna 

### ABSTRACT

This is a protocol of a method used to obtain the total lipid extract (TLE) from archaeological sediment samples and the further elution of the TLE into different fraction and their preparation for GC-MS analysis.

### PROTOCOL STATUS





#### Working

We use this protocol in our group and it is working







#### Preparation

- 1 Oven drying of sediment samples at  60 °C during  48:00:00 and subsampling to  5 g of homogenized sediment

#### Lipid Extraction

- 2 Add  40 ml dichloromethane/methanol (DCM/MeOH) 9:1 v/v to the sediment
- 3 Mix for  00:30:00 in a sonicator (below  30 °C )
- 4 Transfer the solvent into a centrifuge tube and centrifugate for  00:10:00 at 4700 rpm
- 5 Repeat step 2-5 two more times
- 6 Filter centrifuged solvent through annealed glass wool and evaporate with N<sub>2</sub>

#### Prepare silica gel column for SPE extraction

- 7 Add small ball of anniled glass wool to column
- 8 Add  0.1 g of sand (50-70 mesh, previously fired at  450 °C during  10:00:00 ) to column
- 9 Add  1 g of silica (70-230 mesh, previously fired at  450 °C during  10:00:00 ) to column

### Adding lipids to column

- 10 Add lipids reconstituted in DCM with about 9 drops every 10 minutes
- 11 Leave the column covered in aluminium foil for at least 12 hours

### Elution

- 12 Calculate dead volume (DV) using n-hexane
- 13 Elute n-alkanes with 3/8 DV using n-Hexane
- 14 Elute aromatics with 2 DV 8/2 v:v n-hexane/DCM
- 15 Elute ketones with 2 DV DCM
- 16 Elute alcohols with 2DV of 1/1 v:v DCM/Ethyl acetate
- 17 Elute acids and diols with 2 DV of Ethyl acetate
- 18 Elute other compounds with 2 DV of Methanol
- 19 Evaporate all fractions using N<sub>2</sub> at 30 °C

### Preparation of alcohols

- 20 Add 100 µl of N,O-Bis(trimethylsilyl)trifluoroacetamide (BSTFA) + trimethylchlorosilane (TCMS) 99:1 v/v and 100 µl of pyridine to the extract
- 21 Derivatize at 80 °C for 01:00:00
- 22 Dry the sample under N<sub>2</sub>

### Preparation of fatty acids

- 23 Derivatize fatty acids by adding 5 ml of MeOH, 400 µl of H<sub>2</sub>SO<sub>4</sub> to the extract
- 24 Heat at 70 °C for 04:00:00

25 Neutralize with saturated sodium bicarbonate solution

26 Extract three times with  3 ml n-hexane and dry under nitrogen

Adding standard and reconstituting

27 Store at  -20 °C until measurement

28 Reconstitute n-alkanes with 5 $\alpha$ -androstan-3-ol (8mg/L) and  150  $\mu$ l of DCM

29 Reconstitute aromatics with 5 $\alpha$ -androstan-3-ol (8mg/L) and  50  $\mu$ l of DCM

30 Reconstitute ketones with 5 $\alpha$ -androstan-3-ol (8mg/L) and  50  $\mu$ l of DCM

31 Reconstitute alcohols with 5 $\alpha$ -androstan-3-ol (8mg/L) and  40  $\mu$ l of DCM and  10  $\mu$ l of Ethyl acetate

32 Reconstitute acids and diols with Methyl C19:0 (8mg/L) and  40  $\mu$ l of DCM and  10  $\mu$ l of Ethyl acetate

33 Reconstitute other compounds with Methyl C19:0 (8mg/L) and  40  $\mu$ l of DCM and  10  $\mu$ l of Ethyl acetate



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited