



Cell Surface Immunofluorescent Staining of Whole Blood [↗](#)

Kelsey Knight¹

¹BioLegend

Sep 13, 2018 dx.doi.org/10.17504/protocols.io.tjwekpe

Version 3

BioLegend

Working



Kelsey Knight

BioLegend



EXTERNAL LINK

http://www.biolegend.com/media_assets/support_protocol/BioLegend_Surface_Staining_Flow_Protocol_060215.pdf

PROTOCOL STATUS

Working

GUIDELINES






Reagent List:

- Cell Staining Buffer (BioLegend Cat. No. 420201)
- Red Cell Lysis Buffer (BioLegend Cat. No. 420301)
- 7-AAD Viability Staining Solution (BioLegend Cat. No. 420403)
- TruStain FcX™ (anti-CD16/32, BioLegend Cat. No. 101319)
- Human TruStain FcX™ (Fc Receptor Blocking Solution, BioLegend Cat. No. 422301)

References:

Current Protocols in Cytometry (John Wiley & Sons, New York), Unit 6 Phenotypic Analysis.

MATERIALS

NAME	CATALOG #	
 Cell Staining Buffer	420201	by BioLegend
 Red Cell Lysis Buffer	420301	by BioLegend
 7-AAD Viability Staining Solution	420403	by BioLegend
 TruStain FcX™	101319	by BioLegend
 Human TruStain FcX™	422301	by BioLegend

1 Add predetermined optimum concentrations of desired fluorochrome conjugated, biotinylated, or purified primary antibodies to 100 µl of anti-coagulated whole blood.


2 Incubate at room temperature for 15-20 minutes in the dark.

 00:20:00

3 Dilute 10X Red Blood Cell (RBC) Lysis Buffer (BioLegend Cat. No. [420301](#)) to 1X working concentration with DI water. Warm to room temperature prior to use. Add 2ml of 1X RBC lysis solution to whole blood/antibody mixture. Incubate at room temperature for 10 minutes.

 00:10:00

4 Centrifuge at 350 X g for 5 minutes, discard the supernatant.

 00:05:00

- 5 Wash 1X with at least 2 ml of Cell Staining Buffer by centrifugation at 350 x g for 5 minutes.


 00:05:00



REAGENT
Cell Staining Buffer
by [BioLegend](#)
Catalog #: [420201](#)

- 6 If using a purified primary antibody, resuspend pellet in residual buffer and add a previously determined optimum concentration of anti-species immunoglobulin fluorochrome conjugated secondary antibody(e.g. FITC anti-mouse Ig) and incubate in the dark for 15-20 minutes.

 00:20:00

- 7 If using a biotinylated primary antibody, resuspend cell pellet in residual buffer and add a previously determined optimum concentration of fluorochrome conjugated Streptavidin (SAv) reagent (e.g. SAv-PE, BioLegend Cat. No. [405204](#)) and incubate for 15-20 minutes in the dark.  00:20:00

- 8 Repeat step 5.

- 9 Resuspend cells in 0.5 ml Cell Staining Buffer or 0.5 ml 2% paraformaldehyde-PBS fixation buffer.

Tip: For gentler fixation (particularly with tandem fluoros), FluoroFix™ Buffer (Cat. No. [422101](#)) may be used.



REAGENT
Cell Staining Buffer
by [BioLegend](#)
Catalog #: [420201](#)

- 10 Analyze with a Flow Cytometer.

Note: If you are unable to immediately read your samples on a cytometer, keep them shielded from light and in a refrigerator set at 4-8°C. The samples should be resuspended in Cell Staining Buffer. Note that samples should not remain in a fixation buffer for extended periods of time as this can affect fluor conformation and fluorescence.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited