

Cas9 M0386 Mixture

New England Biolabs

Abstract

Citation: New England Biolabs Cas9 M0386 Mixture. protocols.io

dx.doi.org/10.17504/protocols.io.ch3t8m

Published: 07 Oct 2014

Protocol

Step 1.

Nuclease-free water 20 μl

Step 2.

10X Cas9 Nuclease Reaction Buffer 3 μl

Step 3.

300nM sgRNA **3 μl** (30nM final)

Step 4.

1 μM Cas9 Nuclease, S.pyogenes (30nM final)



1 μl Additional info:



Cas9 Nuclease, S. pyogenes - 70 pmol M0386S by New England Biolabs