

Sep 24, 2019

## Marchantia agrobacterium transformation of sporelings in multi-well plates

Linda Silvestri<sup>1</sup>, Eftychis Frangedakis<sup>2</sup>, Susana Sauret-Gueto<sup>1</sup>, Marius Rebmann<sup>1</sup>

<sup>1</sup>Plant Sciences, University of Cambridge, OpenPlant, <sup>2</sup>University of Cambridge, Plant Sciences, OpenPlant

1 Works for me dx.doi.org/10.17504/protocols.io.48tgzwn

OpenPlant Project



Eftychis Frangedakis  
University of Cambridge, Plant Sciences, OpenPlant



### ABSTRACT

A modification of the Ishizaki et al 2008 *Agrobacterium* mediated Marchantia sporeling transformation protocol is used.

Sterilised spores are grown for 5-7 days in ½ strength Gamborg plates prior to co-cultivation for 2 days with agrobacterium in liquid media in multiwell plates. Sporelings are then spread on media with the appropriate selective antibiotic. In about 5 days, positive transformants start to emerge.

#### 1 Day -7 Spore preparation:

Marchantia spores are grown for 7 days in ½ strength Gamborg agar plates (A in Figure).

#### 2 Day -2 Agrobacterium preparation:

- Inoculate a single colony of agrobacterium, transformed with the plasmid of interest, into 5mL LB media plus antibiotics, and then incubate at 28°C with shaking at 150 rpm for two days.

#### Day 1

- After two days, centrifuge the 5mL agrobacterium culture for 10 minutes at 2000g
- Re-suspend in 5ml of ½ strength Gamborg's B5 plus vitamins media with 1% sucrose.
- Add acetosyringone to a final concentration of 100 µM
- Incubate with shaking for 6h at 28°C at 150 rpm.

#### 3 Day 1

Using a sterile scalpel transfer the 5-7 days old Marchantia sporelings in 6 mL of of liquid ½ Gamborg media in a 50 mL Falcon tube (B and C in Figure) and mix well.

#### 4 Aliquote 3 mL of liquid ½ Gamborg B5 vitamins plus supplements media in each well of a 6-well plate.


#### 5 Aliquote 1 mL of spores in each well (D in Figure).

- 6 Add acetosyringone to each well to a final concentration of 100  $\mu$ M.
- 7 Place the 6- well plate on a shaker at 120 rpm for 2 days at 21°C under continuous light (E in Figure).
- 8 Using a pipette transfer the sporeling in a 70 or 100  $\mu$ m cell strainer placed in a 50 mL Falcon tube (F and G in Figure).
- 9 Wash the sporelings with 50 mL of sterile water supplemented with 100  $\mu$ g/ml cefotaxime to remove excessive agrobacterium.
- 10 Plate on  $\frac{1}{2}$  strength Gamborg B5 plus vitamins 1.2% agar plates with cefotaxime and antibiotics (H and I in Figure)

## 11 Day 7

After 5-7 days succesful transformants start to be visible on the plate (J in Figure)

12



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited