

# Immunohistochemistry Protocol for Beta Amyloid Products using USA Detection Kit

Kelsey Miller

## Abstract

**Citation:** Kelsey Miller Immunohistochemistry Protocol for Beta Amyloid Products using USA Detection Kit. **protocols.io**  
dx.doi.org/10.17504/protocols.io.e2ibgce

**Published:** 03 Jun 2016

## Guidelines

-Protocol can be used for Beta Amyloid products that list "IHC" as an application on the datasheet (ex 4G8, 6E10, etc)

-Use with Ultra Streptavidin Detection Kit ([SIG-32250](#)) or ([SIG-32248](#)). All steps should be done in a humidity chamber such as [SIG-31031](#)

## Protocol

### Clear Slides

#### Step 1.

Clear Slides: Removes paraffin and hydrates the tissue

*Note: If using frozen sections, allow slides to come to room temperature for 15 minutes & proceed to step (F) only*

A.	Xylene	5 minutes in each of (3) different 250mL containers
B.	100% alcohol	5 minutes in each of (3) different 250mL containers
C.	95% alcohol	3 minutes in (1) 250mL container
D.	70% alcohol	3 minutes in (1) 250mL container
E.	water	1 minutes in each of (3) different 250mL containers

F. H2O2 (3%) 15 minutes in (1) 250mL container

---

## Rinse Slides

### Step 2.

Rinse slides with lab grade water.

*Note: Lab grade filtered water such as injection grade, cell culture grade, Reverse Osmosis De-ionization (RODI).*

## Antigen Retrieval

### Step 3.

Antigen Retrieval (refer to product datasheet, not always required)

70% Formic Acid – incubate the slides for 20 minutes at room temperature.

*Note: This antigen retrieval step is harsh on the tissue. If using frozen sections reduce time to 5-10 minutes or omit if tissue falls off the slide*

 DURATION

00:20:00

## Antigen Retrieval

### Step 4.

Rinse Slides with 1X PBS

### Step 5.

Apply serum block for at least 5 minutes.

Do NOT wash after this step

 DURATION

00:05:00

### Step 6.

Blot off serum block

### Step 7.

Apply primary antibody – dilute to 1 mg/mL in PBS

#### **Step 8.**

Incubate primary antibody 60 minutes at room temperature

 **DURATION**

01:00:00

#### **Step 9.**

Rinse slides with 1X PBS

#### **Step 10.**

Apply USA Linking reagent - 20 minutes incubation

 **DURATION**

00:20:00

#### **Step 11.**

Rinse slides with 1X PBS

#### **Step 12.**

Apply Labeling Reagent – 20 minutes incubation

 **DURATION**

00:20:00

#### **Step 13.**

Rinse with 1X PBS

#### **Step 14.**

Apply chromogen – 5 minutes incubation.

Dilute according to manufacturer's instructions

1. AEC Chromogen: 20mL AEC chromogen + 1mL AEC substrate buffer
2. DAB Chromogen: 40 mL DAB chromogen + 1mL DAB substrate buffer

 **DURATION**

00:05:00

#### **Step 15.**

Rinse slides with lab grade water

Counterstain

### Step 16.

Submerge slides in Mayer's Hematoxylin for 30 seconds

 DURATION

00:00:30

Counterstain

### Step 17.

Rinse under running lab grade water for 1 minute or until water is clear

 DURATION

00:01:00

Counterstain

### Step 18.

Submerge slides in Bluing Reagent for 1 minute

 DURATION

00:01:00

Counterstain

### Step 19.

Rinse under running lab grade water for 1 minute

 DURATION

00:01:00

Clear slides

### Step 20.

Clear slides: Dehydrate the tissue.

1. 95% alcohol 3 minutes in (1) 250mL container
2. 100% alcohol 5 minutes in each of (3) different 250mL container
3. Xylene 5 minutes in each of (3) different 250mL container

Coverslip

### Step 21.

Cover slip slide using Permanent Aqueous Mounting Medium ([SIG-31010](#)).

*Note: Do not use xylene based mount with AEC Chromogen as it will dissolve the chromogen*