

# Preparation of Proteoliposomes

Zehra Kahveci

## Abstract


Detergent mediated reconstitution of membrane proteins.

**Citation:** Zehra Kahveci Preparation of Proteoliposomes. [protocols.io](https://doi.org/10.17504/protocols.io.rvwd67e)

[dx.doi.org/10.17504/protocols.io.rvwd67e](https://doi.org/10.17504/protocols.io.rvwd67e)

**Published:** 23 Jul 2018

## Materials

 Poly(ethylene glycol) octyl ether 40530 by Sigma Aldrich

 Bio-Beads™ SM-2 Resin 1523920 by BioRad Sciences

## Protocol

### Liposomes preparation

#### Step 1.



. Giant Unilamellar Vesicles (GUVs)

Preparation by Electroformation Method

CONTACT: [Zehra Kahveci](#)

### Preparation of phospholipid stock solution

#### Step 1.1.

Weigh phospholipid of interest, dissolve the powder in chloroform at 10 mM final concentration.

#### Step 1.2.

For fluorescent GUVs, add fluorescent phospholipid Texas-Red-DHPE (1 mM stock) to reach the final molar ratio of 1:500.

### Electroformation by using Nanion Vesicle Prep Pro

#### Step 1.3.

Identify the conductive sides of the ITO-slides with a multimeter.

Electroformation by using Nanion Vesicle Prep Pro

**Step 1.4.**

Use a 5 µl Hamilton syringe to spread 10 µl of phospholipid stock solution on one of the ITO glass slide.

Electroformation by using Nanion Vesicle Prep Pro

**Step 1.5.**

Leave the ITO coating side upwards into a vacuum dessicator for 15 min.

Electroformation by using Nanion Vesicle Prep Pro

**Step 1.6.**

Prepare 195 mM sucrose solution, dissolve it in Hepes Buffer 5 mM pH 7.4, vortex until sucrose is dissolved.

Electroformation by using Nanion Vesicle Prep Pro

**Step 1.7.**

Spread a thin layer of vacuum grease evenly on the O-ring that you use.

Electroformation by using Nanion Vesicle Prep Pro

**Step 1.8.**

After the ITO-slides are dried, take the top part off from the Nanion Vesicle Prep Pro chamber by removing the screws.

Electroformation by using Nanion Vesicle Prep Pro

**Step 1.9.**

Use tweezers to place the ITO-slide lipid film point upwards in the Nanion Vesicle Prep Pro.

Electroformation by using Nanion Vesicle Prep Pro

**Step 1.10.**

Place an O-ring around the dried lipid film, add 280 µl of 195 mM sucrose solution into the O-ring.

Electroformation by using Nanion Vesicle Prep Pro

**Step 1.11.**

Place a second ITO-slide on top of the O-ring with the conductive side facing downwards using tweezers. Note that, the right end of the top slide has to touch properly to the electrode on the right-side.

Electroformation by using Nanion Vesicle Prep Pro

**Step 1.12.**

Place the top part of the chamber back on and screw it carefully.

#### Electroformation by using Nanion Vesicle Prep Pro

##### **Step 1.13.**

Turn-on the Nanion VPP, select the appropriate protocol and initiate the formation process.

##### TEMPERATURE

37 °C Additional info:

##### DURATION

00:05:00 Additional info: Frequency 5 Hz, applied voltage from 0 to 3 V

##### DURATION

02:00:00 Additional info: Frequency 5 Hz, applied voltage 3 V

##### DURATION

00:05:00 Additional info: Frequency 5 Hz, applied voltage from 3V to 0

#### Electroformation by using Nanion Vesicle Prep Pro

##### **Step 1.14.**

After the protocol ends remove the screws and take off the top part of the chamber. Use the tweezers to remove the top ITO-slide.

#### Electroformation by using Nanion Vesicle Prep Pro

##### **Step 1.15.**

Cut a pipette tip and collect the sample from the ITO-slide surface slowly.

#### Fluorecence Microscopy

##### **Step 1.16.**

Withdraw a 5 µl aliquot of the electroformation product and inspect the GUVs under an epifluorescence microscope.

#### Proteoliposomes-DAY1

##### **Step 2.**

Prepare 2 eppendorfs containing 100 µL of 195 mM sucrose solution (freshly made in Hepes buffer 5mM pH 7.4) and add 100 µL of electroformation product into each eppendorf.

#### Proteoliposomes-DAY1

##### **Step 3.**

Shake the Eppendorfs gently upside-down.

#### Proteoliposomes-DAY1

##### **Step 4.**

Dilute the membrane transporter (5mg/mL) in Poly(ethylene glycol) octyl ether (**o-poe**), 1:1 (v:v), total volume of 1 $\mu$ L.

#### Proteoliposomes-DAY1

##### **Step 5.**

Add 1 $\mu$ L of membrane transporter:o-poe 1:1 (v:v) into Eppendorf 1 (**SAMPLE**)

#### Proteoliposomes-DAY1

##### **Step 6.**

Add 1 $\mu$ L of o-poe into Eppendorf 2 (**BLANK**)

#### Proteoliposomes-DAY1

##### **Step 7.**

Shake the **SAMPLE** and **BLANK** gently upside-down.

#### Proteoliposomes-DAY1

##### **Step 8.**

Place both tubes on to a shaker, shake them at 600 rpm.

 **TEMPERATURE**

24 °C Additional info:

 **DURATION**

01:00:00 Additional info:

#### Proteoliposomes-DAY1

##### **Step 9.**

Add 2.5 mg Biobeads into both tubes.

#### Proteoliposomes-DAY1

##### **Step 10.**

Shake the tubes upside-down, place them on to a shaker and shake them at 600 rpm.

 **TEMPERATURE**

24 °C Additional info:

 **DURATION**

01:00:00 Additional info:

#### Proteoliposomes-DAY1

##### **Step 11.**

Store the **SAMPLE** and **BLANK** in the fridge at 4°C, overnight.

#### Proteoliposomes-DAY 2

**Step 12.**

Cut a pipette tip, collect the SAMPLE and the BLANK from each each eppendorf and put the SAMPLE and BLANK into new eppendorfs.

---

**Proteoliposomes-DAY 2****Step 13.**

The SAMPLE and the BLANK can be used in the following 4-5 days. (Store in the fridge).

---