



Jul 27, 2019

## Growing bacteria in b-broth

Cristian Riccio<sup>1</sup>

<sup>1</sup>University of Cambridge

1

Works for me

dx.doi.org/10.17504/protocols.io.5xjg7kn



Cristian Riccio  
University of Cambridge



### GUIDELINES

Work in sterile conditions

### MATERIALS

NAME	CATALOG #	VENDOR
b-broth	<a href="#">View</a>	
Agar plate with single colonies of bacteria	<a href="#">View</a>	
Sterile pipette tips	<a href="#">View</a>	
Bunsen burner	<a href="#">View</a>	
37 C incubator	<a href="#">View</a>	

### SAFETY WARNINGS

Do not wear gloves around the flame of the Bunsen burner.

- 1 Find an isolated colony of bacteria on the agar plate.
- 2 At the Bunsen burner, touch the colony with the tip of a pipette tip and drop the pipette tip into a 100 ml bottle of b-broth.
- 3 Put the 100 ml bottle of b-broth with the bacteria in the incubator, 37 C for *E. coli*.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited