

# Cloning guides to lentiCRISPR v2

Amit Weiner

## Abstract

I have been struggling with cloning guide RNA sequences to the lentiCRISPR vector (<https://www.addgene.org/52961/>) for months now and have finally figured out a method that consistently works for me, so I felt like I should share this information with whomever is struggling as much as I did!

**Citation:** Amit Weiner Cloning guides to lentiCRISPR v2. **protocols.io**

[dx.doi.org/10.17504/protocols.io.qx3dxqn](https://doi.org/10.17504/protocols.io.qx3dxqn)

**Published:** 13 Jun 2018

## Before start

Make Stbl3 chemically competent cells.

## Protocol

### Vector preparation:

#### Step 1.

Digest the lentiCRISPRv2 vector:

3µg vector

2µl 10X Tango buffer (Thermo Fisher)

1µl 20mM DTT

1-1.5µl Esp3I (Thermo Fisher)

Water to 20µl

In a thermocycler:

37°C for 4 hours, inactivate at 65°C for 20 mins, keep at 4°C.

**SKIP** the alkaline phosphatase step.

#### Step 2.

Use a gel purification/PCR cleanup kit (Qiagen) **WITHOUT** running on a gel.

Optional: run 200ng of the purified vector on a gel to verify digestion.

### Guide insert preparation:

#### Step 3.

Anneal and phosphorylate gRNA oligos pair:

1µl of each oligo (100µM stock)

1µl 10X T4 ligation buffer (not PNK buffer)

0.5µl T4 PNK

Water to 10µl

In a thermocycler: 37°C for 30 mins, 95°C for 5 mins, ramp down to 25°C at 0.1°C/sec (or 5-6°C/min).  
Optional: keep at 4°C.

#### Step 4.

Serially dilute the annealed oligos to 1:500

#### Ligation:

#### Step 5.

Out of ligation at a vector:insert molar ratios of 1:5, 1:10, 1:20, I found that 1:5 works best.

50ng vector

1.5µl 10X T4 ligation buffer (NEB)

1µl T4 ligase (NEB)

2µl diluted oligos (1:500)

Water to 15µl

Incubate at RT for 1-2 hours.

#### Transformation:

#### Step 6.

Transform 5µl of the ligation reaction to 50µl Stbl3 chemically competent cells.

If you incubate the Stbl3 cells at 30°C the colonies will be VERY small so look for them carefully. Incubating them at 37°C didn't result in LTR recombination in my hands.

#### Colony PCR:

#### Step 7.

Replica plate colonies on a new LB-Amp plate prior to inserting the tip to 10µl water.

Mix well/ vortex

Boil at 98°C for 10 mins

#### Colony PCR:

#### Step 8.

Using primers upstream and downstream of the guide insert sequence, perform colony PCR:

F primer: gca tat acg ata caa ggc tgt tag aga ga

R primer: gag cca gta cac gac atc act t

	µl
10X buffer	2
10mM dNTPs	0.4
F primer (10uM)	0.3

R primer (10uM)	0.3
Boiled colony	0.5
Taq pol	0.1
ddH <sub>2</sub> O	16.4
Total:	20

\* Can be scaled up or down.

PCR program:

1	94°C	5 min
2	94°C	30 sec
3	54°C	1 min
4	72°C	1 min/Kb
5	Go to step 2	X30
6	72°C	5 min
7	4°C	Hold

#### Colony PCR:

##### Step 9.

Run on a gel -

product of positive colonies: 500 bp

product of negative colonies: 2500 bp

#### Colony PCR:

##### Step 10.

Send 2 colonies per guide for sequencing with U6\_fwd commercial primer.