



May 14,
2019

Working

UC Davis - HbA1c Protocol [↗](#)

Peter Havel¹

¹University of California, Davis

[dx.doi.org/10.17504/protocols.io.ynrfvd6](https://doi.org/10.17504/protocols.io.ynrfvd6)

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org

Lili Liang

ABSTRACT

Summary:

Direct Enzymatic HbA1c test is an enzymatic assay in which lysed whole blood samples are subjected to extensive protease digestion with *Bacillus* sp protease. This process releases amino acids including glycosylated valines from the hemoglobin beta chains. Glycosylated valines then serve as substrates for specific recombinant fructosyl valine oxidase (FVO) enzyme, produced in *E. coli*. The recombinant FVO specifically cleaves N-terminal valines and produces hydrogen peroxide. This, in turn, is measured using a horseradish per-oxidase (POD) catalyzed reaction and a suitable chromagen.

The HbA1c concentration is expressed directly as %HbA1c by use of a suitable calibration curve in which the calibrators have values for each level in %HbA1c.

EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=89&docType=Protocol>

MATERIALS

NAME	CATALOG #	VENDOR	CAS NUMBER	RRID
Calibrator	DZ168A-CAL	Diazyme		
Reagents	DZ168A-K	Diazyme		
Microplate				
Platereader				

MATERIALS TEXT

Reagent Preparation:

Lysis Buffer – ready to use

R1a & R1b – mix together in 70:30 ratio

R2 – ready to use

- 1 Use 250 µl of lysis buffer to lyse 20 µl samples of whole blood and calibrators.

IMPORTANT: Make sure the samples are totally lysed. Any solid material floating around will interfere with reading in the platereader.

- 2 Mix R1a and R1b reagents together in a 70:30 ratio.

3 Add 25 µl of each calibrator and sample to each well.

4 Add 160 µl of reagent R1ab mix to each well. Incubate at 37°C for 5 minutes then read at 720 nm.

IMPORTANT: Make sure not to add any bubbles to the wells when dispensing reagents, this will interfere with reading in the platereader.

5 Add 70 µl of R2 to each well. Incubate at 37°C for 3 minutes then read at 720 nm.

6 Subtract blank readings from final readings. The assay will be linear so the unknown samples can be calculated as $(\text{Sample Absorbance} \div \text{Calibrator Absorbance}) \times \text{Calibrator Concentration}$.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited