

Propagation/Amplification of Lytic Agent

Mathias Middelboe, Amy M. Chan, and Sif K. Bertelsen

Abstract

For use in "[Isolation of cyanophages by liquid enrichment assay](#)"

Citation: Mathias Middelboe, Amy M. Chan, and Sif K. Bertelsen Propagation/Amplification of Lytic Agent. **protocols.io**
dx.doi.org/10.17504/protocols.io.dq85zv

Published: 08 Feb 2016

Protocol

Step 1.

Set up bioassay using 5 mL or larger culture tubes in triplicates.

Step 2.

Add between 5 to 50 μ L of each sample below to target cells in log phase.

📌 NOTES

Amy Chan 03 Sep 2015

- Whole lysate (unfiltered).
- Filtered lysate (0.22 or 0.45 μ m)
- Negative control (no addition, or use filtered media)

Step 3.

Monitor in vivo chlorophyll fluorescence for about 1 week, look for decrease in relative fluorescence compared with control cultures.

Step 4.

If the cultures lyse, then the lytic agent is most probably a virus.

Step 5.

Propagate the lytic agent several times to dilute out nonreplicating viruses.

Step 6.

Filter the lysate and use it to obtain pure clonal stocks.