

# QuikChange Mutagenesis

James Fraser

## Abstract

**Citation:** James Fraser QuikChange Mutagenesis. **protocols.io**

dx.doi.org/10.17504/protocols.io.dij4cm

**Published:** 10 Aug 2015

## Guidelines

### Primers

I use enzyme x to work with DNA sequences, it is installed on all the macs and has nice translation and reverse complement functions. Many of the same functions are available in a plasmid editor , which has Mac, Windows, and Linux versions. Choose a codon for the new amino acid that is the most similar to the original codon.

-Forward Primer

-identify the bases that code for your residue of interest (eg. 295-297)

-A - Copy the 10 bases before the codon (eg. 285-294)

-B - Write your new codon after the 10 bases you just copied

-C - Copy the 22 bases that follow the codon (eg. 298-319), I always terminate in a G or C your forward primer should be 5'-A-B-C-3'

-Reverse Primer

-D - Copy the ten bases after the codon (eg. 298-307) - **then reverse complement it!**

-E - Write your new codon - **then reverse complement it!**

-F - Copy the 22 bases that precede the codon (eg. 273-294) - **then reverse complement it!**  
your reverse primer should be 5'-D-E-F-3'

-Order your primers from ElimBio (The newest PO B000177051)

-They should arrive in 200uM concentration, make a 20uM stock by diluting 5ul of Elim Primer with 45ul of dH2O

### PCR

#### Recipe

-Recipe 1.25ul of 20uM Primer F

-1.25ul of 20uM Primer R

-10ul Phusion Buffer (5x) 1ul of 10mM dNTPs

-0.5ul Template DNA (from miniprep, preferably from DH5alpha cells, but BL21 is still dam+ so it should be fine - 5ng starting template is plenty! (note: good range = 50 - 100 ng) )

-0.5ul Phusion polymerase

-35ul H2O

#### Cycle

-98C for 30s 98C for 5s (I usually do this step for 10 seconds)

-53C for 20s

-72C for 20s/kb (usually plasmids are ~7kb = 2:20)

-Cycle 16 times - *more cycles are actually bad!*

- 72C for 8:00
- 4C for hold

**If it does not work:**

- Miniprep more colonies... ;)
- Order PAGE purified primers if you haven't already. We have a good price from Invitrogen.
- Redesign your primers with netprimer
- For high background of WT sequence, decrease starting amount of template (really, 5ng is plenty),
- Get fresh Dpn1 and digest ~24hrs.
- Run a gradient of annealing temperatures (try 52-62).

Split reaction protocol (for high primer-dimer problems):

**Recipe F**

- 1ul of 20uM Primer F
- 5ul Phusion Buffer (5x)
- 0.5ul of 10mM dNTPs
- 0.5ul Template DNA (from miniprep, preferably from DH5alpha cells, but BL21 is still dam+ so it should be fine) 0.5ul Phusion polymerase
- 17.5 ul H2O

**Recipe R**

- 1ul of 20uM Primer F
- 5ul Phusion Buffer (5x)
- 0.5ul of 10mM dNTPs
- 0.5ul Template DNA (from miniprep, preferably from DH5alpha cells, but BL21 is still dam+ so it should be fine)
- 0.5ul Phusion polymerase
- 17.5 ul H2O

**Cycle each F and R in separate tubes for:**

- 98C for 30s
- 98C for 5s (I usually do this step for 10 seconds)
- 53C for 20s
- 72C for 20s/kb (usually plasmids are ~7kb = 2:20)
- Cycle 10 times
- 72C for 8:00
- 4C for hold

**Mix tubes F and R and cycle:**

- 98C for 30s
- 98C for 5s (I usually do this step for 10 seconds)
- 53C for 20s
- 72C for 20s/kb (usually plasmids are ~7kb = 2:20)
- Cycle 10 times
- 72C for 8:00
- 4C for hold

Follow transformation, DpnI digest as above

**Degenerator:** <http://andersonlab.qb3.berkeley.edu/Software/Degenerator.html>

NB: Not to be confused with this: <http://www.imdb.com/title/tt0100449/>

## Materials

QIAquick PCR Purification Kit [28104](#) by [Qiagen](#)

## Protocol

### PCR

#### Step 1.

Prepare PCR mixture.

#### ✓ PROTOCOL

#### . [PCR Mixture for QuickChange Mutagenesis](#)

CONTACT: [James Fraser](#)

#### Recipe

##### Step 1.1.

125ul of 20uM Primer F.

#### Recipe

##### Step 1.2.

1.25ul of 20uM Primer R.

#### Recipe

##### Step 1.3.

10ul Phusion Buffer (5x).

#### Recipe

##### Step 1.4.

1ul of 10mM dNTPs.

#### Recipe

##### Step 1.5.

0.5ul Template DNA.

#### 🗨 NOTES

**James Fraser** 27 Jul 2015

From miniprep, preferably from DH5alpha cells, but BL21 is still dam+ so it should be fine - 5ng starting template is plenty! (note: good range = 50 - 100 ng).

#### Recipe

##### Step 1.6.

0.5ul Phusion polymerase.

#### Recipe

##### Step 1.7.

35ul H2O.

## PCR

### Step 2.

Perform the PCR using the following cycling conditions:

- 98C for 30s
- 98C for 5s (I usually do this step for 10 seconds)
- 53C for 20s
- 72C for 20s/kb (usually plasmids are 7kb = 2:20)
- Cycle 16 times - more cycles are actually bad!
- 72C for 8:00
- 4C for hold

## DpnI and transformation

### Step 3.

After PCR, add 1ul of DpnI to each pcr tube.

## DpnI and transformation

### Step 4.

Incubate 1 hour-O/N at 37C.

 **DURATION**

15:00:00

 **NOTES**

**James Fraser** 27 Jul 2015

I always do this O/N or for at least 4 hours (Avi).

**James Fraser** 27 Jul 2015

I got 95% efficiency when I DpnI ~24hrs (Cat).

## DpnI and transformation

### Step 5.

PCR Purify using Qaigen kit.

 **REAGENTS**

QIAquick PCR Purification Kit [28104](#) by [Qiagen](#)

## DpnI and transformation

### Step 6.

Add 5ul of PCR reaction into 25ul DH5alpha or TG1 cells.

 **NOTES**

**James Fraser** 27 Jul 2015

I eluted in 27ul EB and transformed all 25ul elution in 50ul TG1 cells + 25ul KCM (Cat).

## DpnI and transformation

### Step 7.

Incubate on ice for 5 minutes.

 **DURATION**

00:05:00

 **NOTES**

**James Fraser** 27 Jul 2015

If you have time, incubate for 30min - 1hr (Avi).

## DpnI and transformation

### Step 8.

Heatshock at 42C for 45 seconds.

 [DURATION](#)

00:00:45

DpnI and transformation

### Step 9.

Recover on ice for 2 minutes.

 [DURATION](#)

00:02:00

DpnI and transformation

### Step 10.

Add 200ul LB (if in 96 well or 8 strip format) to each reaction.

 [NOTES](#)

**James Fraser** 27 Jul 2015

350ul SOC, plate 50ul (Cat).

DpnI and transformation

### Step 11.

Shake at 37C for 1 hour.

 [DURATION](#)

01:00:00

DpnI and transformation

### Step 12.

Plate on warmed antibiotic plates.