

Marine Purity Broth

Bertillson et al.

Abstract

Purpose: To test axenic cultures for purity

Citation: Bertillson et al. Marine Purity Broth. **protocols.io**

dx.doi.org/10.17504/protocols.io.dcm2u5

Published: 02 Feb 2016

Guidelines

Formulation:

17 g BD Difco AC Broth

20 g NaCl

Qs to 1 L with distilled H₂O

After broth has cooled:

2 M MgSO₄ (2 M MgSO₄·7H₂O = 49.3 g/100 ml H₂O, dissolve and 0.2 µm filter)

1 M CaCl₂ (1 M CaCl₂·2H₂O = 14.7 g/100 ml H₂O, dissolve and 0.2 µm filter)

Bertillson et al. 2003 Limnol Oceanogr 48(5): 1721-1731

Protocol

Step 1.

Autoclave BD Difco AC Broth and NaCl Qs to 1 L with distilled water for 20 min.

🕒 **DURATION**

00:20:00

📝 **NOTES**

VERVE Team 12 Aug 2015

See formulation in guidelines.

Step 2.

After broth has cooled, add 16 ml 2 M MgSO₄ and 10 ml 1 M CaCl₂

📝 **NOTES**

VERVE Team 12 Aug 2015

(2 M MgSO₄·7H₂O = 49.3 g/100 ml H₂O, dissolve and 0.2 µm filter)

VERVE Team 12 Aug 2015

(1 M CaCl₂·2H₂O = 14.7 g/100 ml H₂O, dissolve and 0.2 µm filter)

Step 3.

Dispense 2.0 or 4.5 ml per sterile tube

Step 4.

Inoculate sterile MPB with 1:5 (vol:vol) dilution of culture

Step 5.

Incubate at 22°C

📌 NOTES

VERVE Team 10 Jul 2015

Heterotrophs normally grow in less than 2 days. No growth in the broth after 1-2 weeks indicates culture is axenic.