

# Intestine cell dissociation

Andrew Potter

## Abstract

Protocol for human intestine cell dissociation.

**Citation:** Andrew Potter Intestine cell dissociation. **protocols.io**

dx.doi.org/10.17504/protocols.io.nf3dbqn

**Published:** 08 Mar 2018

## Before start

### Checklist prior to beginning:

- Centrifuges, large and small, set to 4 C
- Make enzyme stock; place 2 tubes of enzyme on dry ice.
- Make 0.01% BSA/PBS (50 mL)
- Things you need: petri dishes, clean forceps, razor blade, pipets, 30 µM filters, timer.

## Stock solution for enzyme

- 895 µL DPBS
- 5 µL 0.5 M EDTA (2.5 mM final)

→Add 100 µL enzyme (100 mg/mL) to 900 µL of enzyme stock to make 1X enzyme mix. Add 28 mg of tissue to each 900 µL of enzyme mix.

## Materials

 Please see Guidelines for required materials by Contributed by users

## Protocol

### Step 1.

While excluding as much PBS as possible, **weigh out tissue** using Mettler.

### Step 2.

After weighing out tissue, **transfer to petri dish** on ice and **mince tissue** using grinding motion with razorblade for 2-3 minutes.

### Step 3.

After tissue is minced finely, **add 1 mL enzyme mix per 28 mg of tissue** to the petri dish and pipet minced tissue + enzyme into eppendorf tube (on ice).

#### AMOUNT

1 ml Additional info: enzyme mix per 28 mg of tissue

### Step 4.

**Start timer.** Leave tube on ice - **initially shake vigorously to break up the tissue**, 3-

5x every 30-45 seconds for 5 minutes.

### Step 5.

Now, when big chunks are broken up, shake every 1 minute while leaving on ice for **10 minutes**.

### Step 6.

**Triturate 10X** with 1 mL pipet set to 700  $\mu$ L.

### Step 7.

Continue to shake vigorously every minute to re-suspend tissue for **15 minutes**.

### Step 8.

**Triturate digest mix again 10X** and spin digest mix at 90 G for 30 seconds.

### Step 9.

Remove supernatant (80%) containing single cells and filter using 30  $\mu$ M filter while leaving chunks on bottom; rinse filter with 10 mL PBS/BSA.

#### AMOUNT

10 ml Additional info: PBS/BSA

### Step 10.

To residual chunks of tissue add additional 1 mL of enzyme (per 28 mg tissue).

#### AMOUNT

1 ml Additional info: enzyme (per 28 mg tissue)

### Step 11.

Shake vigorously 3-4X every minute for **10 additional minutes**.

### Step 12.

**Triturate again 10X.**

### Step 13.

Continue to shake vigorously every minute for **15 minutes** total additional time.

### Step 14.

Triturate again 10X and filter using the same 30  $\mu$ M filter and rinse with 10 mL PBS/BSA.

#### AMOUNT

10 ml Additional info: PBS/BSA

### Step 15.

Divide flow-through into 2 15 mL tubes.

### Step 16.

Spin 600 g for 5 minutes at 4 °C.

#### TEMPERATURE

4 °C Additional info: Spinning

### Step 17.

Carefully remove and save supernatant - re-suspend in 1 mL total PBS/BSA in 1.5 mL tube.

#### AMOUNT

1.5 ml Additional info: PBS/BSA

### Step 18.

Spin down 600 G for 5 minutes - remove and save supernatant.

### Step 19.

**Re-suspend** in 700  $\mu$ L RBC lysis buffer +100  $\mu$ L PBS/BSA (800  $\mu$ L total).

#### AMOUNT

700  $\mu$ L Additional info: RBC lysis buffer

#### AMOUNT

100  $\mu$ L Additional info: PBS/BSA

### Step 20.

Incubate for 3 minutes on ice.

### Step 21.

Spin 600 G for 5 minutes.

### Step 22.

Remove supernatant.

### Step 23.

Briefly re-suspend and check cells in a small volume of PBS/BSA to ensure that there are no more RBCs present.

### Step 24.

**Re-suspend** in 1.3 mL total PBS/BSA.



1.3 ml Additional info: PBS/BSA

### Step 25.

Spin 600 g for 5 minutes at 4 °C.



4 °C Additional info: Spinning

### Step 26.

Remove supernatant and re-suspend in a small volume to check cell concentration.

### Step 27.

Analyze quantity and viability of cells using a hemocytometer with trypanblue; **for Chromium, make concentration to 1 million cells per mL.**