

Genomic DNA extraction from cells

Bao Thai

Abstract

Citation: Bao Thai Genomic DNA extraction from cells. [protocols.io](https://doi.org/10.17504/protocols.io.nk4dcyw)

[dx.doi.org/10.17504/protocols.io.nk4dcyw](https://doi.org/10.17504/protocols.io.nk4dcyw)

Published: 03 Mar 2018

Protocol

Step 1.

Extract genomic DNA from cells with Quick Extract buffer (QE - Epibio #QE09050).

Step 2.

Remove media.

Step 3.

Trysinize cells.

Step 4.

Centrifuge at 16K x g for 2 minutes.

Step 5.

Wash with ice cold 1xPBS.

Step 6.

Centrifuge at 16K x g for 2 minutes.

Step 7.

Add 100ul QE / well (24-well plates, scale up for your plate size).

Step 8.

Let sit 5 minutes.

Step 9.

Transfer to PCR tubes, 100 uL each tube.

Step 10.

Vortex for 15sec, centrifuge to collect liquid in bottom.

Step 11.

Put PCR tubes in thermocycler: 20min at 65C, 20min at 95C, store at -20C.

Step 12.

Measure concentration with Nanodrop. Use QE buffer as blank.