

RNA extration for the Betta splendens genome

Xin Liu

Abstract

This protocol is used to clarify the process of RNA extration for our Betta splendens genome.

Citation: Xin Liu RNA extration for the Betta splendens genome. **protocols.io**

dx.doi.org/10.17504/protocols.io.qq7dvzn

Published: 08 Jun 2018

Protocol

Sample preparation

Step 1.

1. Pour 1.5ml TRIZOL reagent.
2. For tissue samples, grind about 60mg with liquid nitrogen into powder and transfer the powder samples into the 2 ml tube contain of 1.5ml Trizol reagent.



REAGENTS

TRIZOL reagent 15596-026 by [Invitrogen - Thermo Fisher](#)

Tissues lysis

Step 2.

Homogenize 2 minutes and place the sample at rest horizontally for 5 minutes to permit the complete dissociation of nucleoprotein complexes.

Phase separation

Step 3.

1. Centrifuge at 12000×g for 5 minutes at 4°C. Transfer the supernatant to a new 2.0ml tube, add 0.3 ml of Chloroform / isoamyl alcohol(24:1) per 1.5 ml of Trizol reagent. Shake the tubes vigorously for 15 seconds.
2. Centrifuge at 12000×g for 10 minutes at 4°C. After centrifugation, the mixture should separates into three layers: the lower phenol-chloroform phase, an interphase, and an upper aqueous phase. RNA remains in the aqueous phase.



REAGENTS

Chloroform / isoamyl alcohol(24:1) 319988/W205702 by [Sigma](#)

RNA precipitation

Step 4.

1. Transfer the aqueous phase to a new 1.5mL tube; add equal volume of supernatant of isopropyl alcohol. Mixing well and place at -20°C for 2 hours for precipitation.
2. Centrifuge at 13600rpm for 20 minutes at 4°C and remove the supernatant.



REAGENTS

isopropyl alcohol W292907 by [Sigma](#)

RNA washing

Step 5.

1. Wash the RNA pellet with 1 ml 75% ethanol. Re-suspend the pellet and centrifuge at 13600rpm for 3 minutes at 4°C. Repeat this step again, Completely remove the ethanol without disturbing the pellet.
2. Air-dry the RNA pellet in the biosafety cabinet.

Dissolve RNA

Step 6.

Add 50µL of DEPC-treated water to dissolve the RNA pellet.



REAGENTS

DEPC-treated water AM9915G by [Ambion](#)