

Methyltransferase Assay

David Dunigan and Irina Agarkova

Abstract

Citation: David Dunigan and Irina Agarkova Methyltransferase Assay. **protocols.io**

dx.doi.org/10.17504/protocols.io.etjbekn

Published: 13 Jun 2016

Guidelines

Reaction Conditions and Buffer, Materials:

- Reaction conditions: 50 mM Tris-HCl, pH 7.8, 10 mM EDTA, 0.5 mM DTT, 2 μ Ci Adenosyl-L-Methionine, S-Methyl (3 H), 3 μ g DNA
- Preparation of reaction buffer, 10X
 - 3000.0 μ L of 1 M Tris-HCl, pH 7.8
 - 1200.0 μ L of 500 mM EDTA, pH 8.0
 - 1500.0 μ L of d-H₂OTo 95.0 μ L of the reaction buffer, add 5.0 μ L of 100 mM DTT.
Use 5.0 μ L per 50.0 μ L reaction.
- Proteinase K, 2 mg/mL (autodigested at 37°C for 60 min)
- ss salmon sperm DNA, 5 mg/mL
- 20% TCA
- 5% TCA

- 500 mM NaOH

Protocol

Step 1.

Set up the assays, to contain 50.0 μL total per assay with 3 μg of DNA and 2 μCi of Adenosyl-L-Methionine, S-Methyl (^3H), and the enzyme to be tested.

Step 2.

Incubate the assays at 37°C for 2 hours.

 DURATION

02:00:00

Step 3.

Add 30 μL (60 μg) of proteinase K to each assay.

Step 4.

Incubate at 65°C for 20 min.

 DURATION

00:20:00

Step 5.

Add 4 μL (20 μg) of ss salmon sperm DNA as a carrier.

Step 6.

Add an equal volume (84 μL) of 20% TCA and incubate at 4°C for 10 min to precipitate the samples.

 DURATION

00:10:00

Step 7.

Centrifuge the assays in the microfuge for 5 min at 4°C.

 DURATION

00:05:00

Step 8.

Discard the supernatants.

Step 9.

Wash the pellets 2X with 200 μL of 5% TCA in the microfuge for 5 min at 4°C.

 DURATION

00:05:00

Step 10.

Resuspend the pellets with 500 μL each of 500 mM NaOH.

Step 11.

Incubate the assays at 65°C for 20 min.

 DURATION

00:20:00

Step 12.

Precipitate the samples with 500 μL of 20% TCA for 10 min at 4°C.

 DURATION

00:10:00

Step 13.

Centrifuge the samples in the microfuge for 5 min at 4°C.

 DURATION

00:05:00

Step 14.

Wash the pellets 2X with 200 µL of 5% TCA in the microfuge for 5 min at 4°C.

 DURATION

00:05:00

Step 15.

Resuspend the pellet with 500 µL of 5% TCA.

Step 16.

Incubate in boiling water for 30 min (loosen the caps or puncture the caps of the tubes).

 DURATION

00:30:00

Step 17.

Centrifuge the samples in the microfuge for 5 min at 4°C to remove the insoluble material.

 DURATION

00:05:00

Step 18.

Count the supernatants in 10 mL of triton-toluene counting solution.