

# Immunoprecipitation assays(about S100A6 in SW480 Cells)

Shanshan Feng

## Abstract

**Citation:** Shanshan Feng Immunoprecipitation assays(about S100A6 in SW480 Cells). **protocols.io**  
dx.doi.org/10.17504/protocols.io.mqdc5s6

**Published:** 18 Jan 2018

## Protocol

### Cell Lysis

#### Step 1.

1. Rinse a 60 mm culture dish of confluent cells with PBS
2. Lyse the cells with 0.5 ml cold Immunoprecipitation Buffer containing 1% SDS, 1 mmol/L of Na<sub>3</sub>VO<sub>4</sub>, 0.1 mol/L of Tris (pH 7.4), and protease inhibitors (10 mg/L of leupeptin, 5 mg/L of aprotinin, 20 mg/L of soybean trypsin inhibitor, and 1 mmol/L of phenylmethylsulfonyl fluoride).
3. Maintain constant agitation for 30 minutes at 4°C.
4. Scrape the cells from the dish.
5. Sonicate on ice for 5 seconds; repeat 4 times
6. Centrifuge for 5 minutes at 4°C
7. Assay for total protein then adjust concentration to approximately 1 mg/ml with Immunoprecipitation Buffer

### Immunoprecipitation

#### Step 2.

8. In a 1.5 ml microcentrifuge tube, add 20 µL of protein A+G Agarose Beads and transferred to a fresh 1.5 mL tube.
9. centrifuging at 1500 rpm for 30 sec at 4 °C and washed three times with 500 µL of lysis buffer.
10. Carefully pipette to remove supernatant
11. Add 2 µg of antibody or 2 mg of IgG to crude cell lysate
11. Incubate overnight at 4 °C.
12. Wash with 500 µl of Immunoprecipitation Buffer by gentle vortex and remove supernatant and discard repeated three times

13. Resuspend bead pellet in 20  $\mu$ l of 1X SDS Sample Loading Buffer
14. Incubate sample at 70°C for 5 minutes.
15. analysis on SDS-PAGE gel and electrophoresis.