



Oct 24, 2019

Electrode cleaning solution

Michael Economo¹¹Boston University

1

Works for me

dx.doi.org/10.17504/protocols.io.8pwhvpe



Michael Economo ⚡

ABSTRACT

How to make a cleaning solution for cleaning silicon probes

MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Sodium chloride	S3014	Sigma Aldrich
Sodium dodecyl sulfate solution	71736-100ML	Sigma Aldrich
0.2 M Tris Buffer pH 8.5	T1503-100G	Sigma Aldrich
Trypsin solution from porcine pancreas	T4549-100ML	Sigma Aldrich

- 1 Make 1.2X strength cleaning solution (167.67 mL)

First, measure **116.67 ml** H₂O in beaker

Add:

0.117 g NaCl**10 ml** 10% SDS**50 ml** 0.2 M Tris Buffer, pH 8.5

- 2 [On day of use] Add **2 ml** 10X trypsin to **10 ml** 1.2X cleaning solution to make 12 mL 1x cleaning solution



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited