



Feb 25, 2019

Working

## Tissue lysis and digestion for MS analysis

David Kotol<sup>1</sup><sup>1</sup>Department of Protein Science, Royal Institute of Technology, Stockholm, Sweden[dx.doi.org/10.17504/protocols.io.yguftww](https://doi.org/10.17504/protocols.io.yguftww)

Human Protein Atlas

David Kotol  
Science fo Life Laboratory / KTH - Royal Institute of Techno...

### ABSTRACT

Here, a workflow for sample preparation and quantification of brain and pancreatic tissue proteins with use of heavy labelled protein standards (QPrESTs) is described. The standards are added before the digestion ensuring high reproducibility and elimination of nuances appearing during proteolytical cleavage of proteins and further steps such as solid phase extraction. QPrESTs are heavy labelled recombinant protein fragments covering regions of proteins with low homology and span over 50-150 aminoacids. They were produced as a part of Human Protein Atlas project in auxotrophic E.coli and individually purified, serving as an excelent resource for MS-based protein quantification.

### PROTOCOL STATUS

#### Working

We use this protocol in our group and it is working

### MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Methanol (MeOH)	10499560	Fisher Scientific
Acetonitrile (ACN)	10660131	Fischer Scientific
Sodium deoxycholate (SDC)	30970	Sigma Aldrich
Urea	U5378	Sigma Aldrich
Thiourea	T8656	Sigma Aldrich
DL-Dithiothreitol (DTT)	43815	Sigma Aldrich
2-Chloroacetamide (CAA)	22790	Sigma Aldrich
EDTA-free Protease Inhibitor Cocktail	11836170001	Sigma Aldrich
Acetone	34850	Sigma Aldrich
Pierce Trypsin Protease	90057	Thermo Fisher Scientific
Formic acid (FA)	15657520	Sigma Aldrich
Trifluoroacetic acid (TFA)	74564	Sigma Aldrich
Triethylammonium bicarbonate (TEAB)	T7408	Sigma Aldrich
Bio-Rad Protein Assay	5000001	BIO-RAD
Ammonium hydroxide	L13168	Alfa Aesar

### STEPS MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Bio-Rad Protein Assay	5000001	BIO-RAD

## MATERIALS TEXT

Buffers were prepared using MilliQ water if not stated differently

### Tissue lysis and digestion

#### Lysis buffer, 1 mL (7M Urea, 2M Thiourea, 2% SDC)

m (Urea) = 0.42 g

m (Thiourea) = 0.15 g

m (SDC) = 0.02 g

1/5 tablet of protease inhibitors

#### 300mM DTT, 500 µL

m (DTT) = 23 mg

#### 500mM CAA, 1 mL

m (CAA) = 46.75 mg

Dissolved in 100mM TEAB

#### 100mM TEAB, 1 mL

100 µL 1M TEAB

#### Solvent A, 100 mL (3% ACN, 0.1% FA)

3 mL 100% ACN

100 µL 100% FA

### Solid Phase Extraction using Strong Cation Exchange

#### Wash Buffer, 1 mL (30% MeOH, 0.1% FA)

300 µL 100% MeOH

1 µL 100% FA

#### Elution Buffer, 1 mL (30% MeOH, 1.25% NH<sub>4</sub>OH)

300 µL 100% MeOH

45 µL 28% NH<sub>4</sub>OH

## SAFETY WARNINGS

Chloroacetamide (2-chloroacetamide) is a chlorinated organic compound used for alkylating reduced cysteine residues and is suspected of damaging fertility.

### Protein extraction

- 1 Incubate adapter of the Tissue Lyser on dry ice for ⌚ 00:30:00



Tissue Lyser LT

QIAGEN 85600 [↗](#)

Bead mill

- 2 Add one 3mm bead to the tissue sample and incubate on dry ice for ⌚ 00:30:00



Tungsten Carbide Beads, 3 mm

QIAGEN 69997 [↗](#)

Beads for use with Tissue Lyser systems

3 Disrupt the tissue using the Tissue Lyser for [⌚ 00:02:00](#)



Tissue Lyser LT

QIAGEN 85600 [↗](#)

Bead mill

4 Add [🧴 250 µl](#) of Lysis Buffer (7M Urea, 2M Thiourea, 2% SDC)

5 Homogenize the tissue using the Tissue Lyser for [⌚ 00:02:00](#)



Tissue Lyser LT

QIAGEN 85600 [↗](#)

Bead mill

6 Centrifuge for [⌚ 00:30:00](#) at 20,000 × g at [🌡 4 °C](#)


7 Precipitate the proteins with 4 volumes of ice-cold acetone and incubate in [🌡 -20 °C](#) over night


8 Centrifuge at 20,000 × g at [🌡 4 °C](#) for [⌚ 00:30:00](#) and discard the supernatant

9 Wash the pellet two times with [🧴 400 µl](#) of ice-cold acetone followed by centrifugation at 20,000 × g at [🌡 4 °C](#) for

10 Remove the acetone excess and dry the pellet

11 Dissolve the pellet in  200  $\mu$ l of Pellet Solubilizer (7M Urea, 2M Thiourea)

12 Dilute 8 times with  100 Milimolar (mM) TEAB and measure protein concentration using Bio-Rad Protein assay according to the manufacturers instructions



**Bio-Rad Protein Assay**  
 by [BIO-RAD](#)  
 Catalog #: 5000001

#### Protein Digestion

13 Take volume that corresponds to  15  $\mu$ g of proteins and add QPrESTs mixture

14 Add DTT to a final concentration of  10 Milimolar (mM) and incubate 🕒 01:00:00 at 🌡 30 °C

15 Add CAA to a final concentration of  50 Milimolar (mM) and incubate 🕒 00:30:00 in the dark at room temperature

16 Add  300 ng of Pierce Trypsin protease and incubate over night at 🌡 37 °C

17 Quench the digestion with TFA to a final concentration of 0.5% (v/v)

#### Solid phase extraction using Strong Cation Exchange

18 Prepare one StageTip per sample by inserting 3 layers of Empore Cation Extraction Discs in a 200uL tip



Empore Discs  
47mm cation extraction disks  
3M 2251-1 [↗](#)

- 19 Activate membrane with  50 µl 100% MeOH, spin at 1000 x g for  00:01:00
- 20 Equilibrate StageTip with  50 µl Wash Buffer (30% MeOH, 0.1% FA), spin at 1000 x g for  00:01:00
- 21 Apply sample to the StageTip, spin at 1000 x g for  00:01:00
- 22 Wash two times with  30 µl Wash Buffer (30% MeOH, 0.1% FA), spin at 1000 x g for  00:01:00  00:01:00
- 23 Elute two times with  20 µl Elution Buffer (30% MeOH, 1.65% NH<sub>4</sub>OH), spin at 1000 x g for  00:01:00  00:01:00
- 24 Vacuum dry for  00:30:00 at  42 °C and store at  -20 °C until MS-analysis
- 25 For MS analysis, dissolve samples in Solvent A and inject 1.75 µg



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited