

# Western Blot

**Angela Tsai**

## Abstract

To identify specific amino acid sequences of a protein, or tag, using fluorescentlytagged antibodies

**Citation:** Angela Tsai Western Blot. **protocols.io**

[dx.doi.org/10.17504/protocols.io.nvxde7n](https://doi.org/10.17504/protocols.io.nvxde7n)

**Published:** 17 Mar 2018

## Protocol

### Step 1.

Assemble gel polymerizing station according to manufacturer's protocol.

### Step 2.

To create a running gel combine:

- a. 2.4 mL dH<sub>2</sub>O
- b. 1.5 mL 1.5M Tris HCl pH 8.8
- c. 2 mL 30% Acrylamide/Bis solution (Bio-Rad)
- d. 60 µL 10% SDS
- e. 60 µL 10% APS (must be made fresh every time)
- f. 6 µL TEMED

### Step 3.

Pour into gel cassette and let polymerize.

### Step 4.

To create a stacking gel combine:

- a. 6.85 mL dH<sub>2</sub>O
- b. 2.5 mL 0.5M Tris HCl pH6.8
- c. 1.3 mL 30% Acrylamide/Bis Solution (Bio-Rad)
- d. 100 µL 10% SDS

e. 100  $\mu$ L 10% APS

f. 10  $\mu$ L TEMED2

**Step 5.**

Pour stacking gel on top of polymerized running gel into gel cassette. Insert gel comb into stacking gel and let polymerize.

**Step 6.**

After polymerization, remove gel comb to generate wells and transfer gel cassettes to gel running dock as per manufacturer's instructions. Fill empty compartments with 1X Running Buffer.

**Step 7.**

Dilute desired protein samples 1:1 with 2X SDS running buffer and load into gel. Add ladder if desired.

**Step 8.**

Allow to run at 140V for 1:30 hours or until dye front is no longer visible.

**Step 9.**

Fill a basin with methanol and a second basin with 1X Trans-Blot Turbo Transfer Buffer.

**Step 10.**

Let PVDF membrane sit in methanol for 2 minutes, then transfer to transfer buffer basin with foam pads. Let sit for 2 minutes.

**Step 11.**

Assemble transfer according to manufacturer's protocol.

**Step 12.**

Transfer proteins from polyacrylamide gel to PVDF membrane using Turbo-Blot Transfer System.

**Step 13.**

Block 1 hour in 5% Milk

**Step 14.**

Incubate overnight with primary antibody in 5% milk

**Step 15.**

Wash 3 times with TBS-Tween

**Step 16.**

Incubate for 1 hour with secondary antibody in 5% milk

**Step 17.**

Wash 3 times with TBS-Tween

**Step 18.**

Develop the blot by preparing the ECL developing solution. Add 1mL of ECL solution to the membrane and allow to incubate for 2 minutes on each side.

**Step 19.**

Image using the Bio-Rad ChemiDoc Gel Imaging System