



Oct 16, 2019

Promega Glucose-Glo Assay

NUS iGEM¹

¹National University of Singapore

1

Works for me

dx.doi.org/10.17504/protocols.io.8aehsbe



NUS iGEM

National University of Singapore



MATERIALS

NAME ▾

CATALOG # ▾

VENDOR ▾

Glucose-Glo(TM) Assay, 50ml

J6022

Promega

1 Prepare glucose detection reagent mix by adding:

1. Luciferin Detection Solution
2. Reductase
3. Reductase Substrate
4. Glucose Dehydrogenase
5. NAD



These reagents are part of a boxed set - Promega's Glucose-Glo Assay.

2 Dilute **2 µl** overnight culture in **198 µl** 1x PBS and mix

3 Aliquot **50 µl** of each sample into three wells of 96-well plate to create triplicates

4 Add **50 µl** of glucose detection reagent mix into each well

5 Include PBS-only wells as blank



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited