



Jan 24,
2020

High-quality RNA purification with on-column DNase treatment from tissue specimens

Magda Bletsa¹, Antonios Fikatas¹, Yiqiao Li¹, Sophie Gryseels^{1,2}, Jan Felix Drexler³, Philippe Lemey¹

¹Department of Microbiology, Immunology and Transplantation, Rega Institute, KU Leuven – University of Leuven, Leuven, Belgium, ²Department of Ecology and Evolutionary Biology, University of Arizona, Tucson, AZ, USA, ³Charité-Universitätsmedizin Berlin, corporate member of Freie Universität Berlin Humboldt-Universität zu Berlin and Berlin Institute of Health, Institute of Virology, Berlin, Germany

1 Works for me [dx.doi.org/10.17504/protocols.io.8ufhwtm](https://doi.org/10.17504/protocols.io.8ufhwtm)



Magda Bletsa
Rega Institute, KU Leuven



ABSTRACT

This protocol can be used for total RNA purification from tissue samples. It includes an intermediate on-column DNase treatment, which results in higher RNA yield and purity. This method is mainly optimised to obtain hepatitis C virus RNA extracts for whole-genome sequencing, but it can also be used for any viral RNA purification from samples with lower viral loads.

GUIDELINES

- During tissue handling, all procedures should be carried out as quickly as possible.
- Only RNA molecules > 200 nucleotides are purified.
- Do not overload the RNeasy spin column (maximum capacity of **700 µl**).
- Buffer RLT may form a precipitate upon storage. Re-dissolve by putting the bottle under warm water for a few minutes.
- Perform all steps at **Room temperature**, unless otherwise stated.
- Always use new collection tubes to eliminate any possible contamination.

MATERIALS

NAME	CATALOG #	VENDOR
RNeasy Mini Kit	74104	Qiagen
RNase-Free DNase Set	79254	Qiagen
Precellys CK28 Lysing Kit Hard Tissue Homogenizing Bertin	10144-516	VWR International

STEPS MATERIALS

NAME	CATALOG #	VENDOR
Precellys CK28 Lysing Kit Hard Tissue Homogenizing Bertin	10144-516	VWR International

BEFORE STARTING

- When using Buffer RPE for the first time, add 4 volumes of ethanol (100%).
- Prepare DNase I stock solution of RNase-free DNase set as described:

Do not open the glass vial. Transfer **550 µl** of the RNase-free water provided into a **1.5 ml** Eppendorf tube. Using a needle and a syringe inject the water into the lyophilised DNase I glass vial. Mix gently by inverting the vial and make sure the powder on the sides of the vial is all well dissolved. **DO NOT VORTEX**.






For long-term storage of DNase I, remove the stock solution from the glass vial, using the syringe. Divide it into **100 µl** aliquots and store at **-20 °C** for up to 9 months. Thawed aliquots can be stored at **4 °C** for up to 6 weeks. Do not refreeze the aliquots after thawing.

Sample homogenisation

- 1 Add  **600 µl** RLT Buffer to the Precellys lysate tubes.






Precellys CK28 Lysing Kit Hard Tissue
Homogenizing Bertin
by VWR International
Catalog #: [10144-516](#)


- 2 Excise a lentil-sized piece of tissue (maximum amount of  **20 mg** for RNAlater stabilized tissues and  **30 mg** for fresh or frozen tissues) and transfer it quickly to the lysate tubes. Make sure that all tissues are immersed into the RLT reagent.
- 3 Place tubes on dry ice for  **00:02:00** (or until frozen) and then thaw quickly.
- 4 Immediately after thawing disrupt the tissues using a conventional rotor-stator homogeniser (Minilys). The recommended lysis should be performed at medium speed ( **4000 rpm**) for  **00:02:00** .



Minilys Personal Homogeniser
Tissue homogeniser
Bertin Instruments P000673-MLYS0-A [↗](#)

- 5 Place tubes on dry ice for  **00:02:00** (or until frozen) and then thaw quickly.
Check the results and repeat the homogenisation until no more visible fragments are present.
- 6 Centrifuge the lysate for  **00:03:00** at full speed. Carefully remove the supernatant by pipetting, and transfer it to a new  **1.5 ml** Eppendorf tube.

RNA binding

- 7 Add  **600 µl** 70% ethanol to the supernatant, and mix immediately by pipetting 5 times. **DO NOT VORTEX OR CENTRIFUGE**. Proceed immediately to the next step.

- 8 Transfer **700 µl** of the sample to an RNeasy spin column placed in a **2 ml** collection tube. Centrifuge for **00:00:30** at **10000 rpm**.
- If the sample is more than **700 µl**, change collection tubes and transfer the rest to the spin column and centrifuge again. Discard the collection tubes and replace with new ones.
- 9 Add **350 µl** Buffer RW1 to the RNeasy spin column and centrifuge for **00:00:30** at **10000 rpm** to wash the spin column membrane. Change collection tubes.

DNase treatment

- 10 Prepare the DNase treatment master mix, as described:
For each sample, add **10 µl** DNase I stock solution to **70 µl** Buffer RDD. Mix by gently inverting the tube, and spin down briefly. **DO NOT VORTEX**.
- 11 Add the **80 µl** DNase I incubation mix directly to the RNeasy spin column membrane, and place on benchtop at **Room temperature** for **00:15:00**.

Washing steps





- 12 Add **350 µl** Buffer RW1 to the RNeasy spin column and centrifuge for **00:00:30** at **10000 rpm** to wash the spin column membrane. Change collection tubes.
- 13 Add **500 µl** Buffer RPE to the RNeasy spin column and centrifuge for **00:00:30** at **10000 rpm**. Change collection tubes.
- 14 Add **500 µl** Buffer RPE to the RNeasy spin column and centrifuge for **00:02:00** at **10000 rpm**. Change collection tubes.

Dry centrifugation

- 15 Place the RNeasy spin column in a new collection tube and centrifuge at full speed for **00:05:00**. If there is still much liquid passing through the column, change collection tubes and centrifuge for **00:01:00**.

Elution and storage

- 16 Place the spin column in a new **1.5 ml** Eppendorf tube. Add **50 µl** RNase-free water directly to the spin column membrane. Incubate for **00:05:00** at **Room temperature**. Then centrifuge for **00:01:00** at **10000 rpm** to elute the RNA.

- 17 Repeat step 16. Use a new  **1.5 ml** tube for another elution round of  **50 µl** .
- 18 Store at --  **20 °C** . Optionally, you can put the column back to the corresponding collection tube and store them at --  **20 °C** for another future elution. Even after the 3rd or 4th elution time, there are still some RNA molecules passing through the column.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited