



Aug 13,  
2019

## RNA Isolation

Natalie Brown<sup>1</sup>

<sup>1</sup>Joslin

**1** Works for me

[dx.doi.org/10.17504/protocols.io.6hjhb4n](https://doi.org/10.17504/protocols.io.6hjhb4n)

**Natalie Brown**  
Joslin Diabetes Center

### ABSTRACT

#### Real Time Polymerase Chain Reaction

RT PCR.docx

- 1 ● Get plates, aspirate media
- 2 ● Add 2mL ice cold DPBS
- 3 ○ Aspirate DPBS
- 4 ● Add 500uL Trizol/well in chemical hood
- 5 ● Bring yellow tubes to hood, collect trizol and add to tubes
- 6 ● Add 100uL chloroform to tubes
- 7 ○ hand shake 15s
- 8 ○ Place on ice 5min until see layers (aq/lipids/bottom)
- 9 ● Spin full speed 20min cold centrifuge
- 10 ○ During, make 70% ethanol on ice (3mL RNase-free water + 7mL ethanol)

- 11 ○ Make new ep set 1-9
- 12 ● In hood, collect top layer after spinning and place into eppendorfs
- 13 ○ \*\*take note of first volume and make all the same
- 14 ○ Mix 1:1 ethanol with supernatant, use same tip to transfer from eppendorfs into pink spin tubes
- 15 ● Spin 9000x 1min bench centrifuge
- 16 ○ Discard flowthrough in hood, keep tube
- 17 ● Add 350uL RW1 buffer to tubes
- 18 ○ Spin 1 min
- 19 ○ During, prepare DNase soln (with syringe)
- 20 ● Add 70uL DNase soln to middle of each column, leave 15min RT
- 21 ● Repeat RW1 and spin
- 22 ○ Discard collection tube and get new
- 23 ● Wash RPE 500uL
- 24 ○ Spin 1min
- 25 ○ Discard collection and keep tubes
- 26 ● Repeat RPE, spin 2min

- 27 ○ Discard tubes
- 28 ● Spin 1min empty to dry
- 29 ○ Get new eppendorfs, label
- 30 ● Place pink tubes on top of eppendorfs
- 31 ● Add 30uL RNase free water, centrifuge 1min
- 32 ○ Close tubes as quickly as possible and keep on ice→ nanodrop



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited