



Working

Extraction of alarm substance in zebrafish

Rhayra Xavier do Carmo Silva¹, Sueslene Prado Rocha², Monica Gomes Lima-Maximino², [Caio Maximino](#)³

¹Universidade Federal do Pará, ²Universidade do Estado do Pará, ³Universidade Federal do Sul e Sudeste do Pará

Sep 20, 2018 dx.doi.org/10.17504/protocols.io.tr3em8n

[Fish behavior and physiology](#)



Caio Maximino

Universidade Federal do Sul e ...



ABSTRACT

Alarm substance, or *Schreckstoff*, is a mixture of compounds produced by specialized club cells in the epidermis of Ostariophysan fish. When lesioned, these cells release the substance, which acts as an aversive olfactory stimulus capable of evoking intense behavioral responses (alarm reaction) in other fish in the shoal (Maximino et al. 2018). Alarm substance is used as an aversive stimulus to study fear-like behavior in Ostariophysan - especially in zebrafish. The present protocol describes a technique to easily extract alarm substance for acute exposure.

TAGS

Zebrafish

Alarm substance

Show tags

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Maximino, C., Silva, R. X. do C., Campos, K. dos S., Oliveira, J. S. de, Rocha, S. P., Pyterson, M. P., ... Maximino, M. L. (2018). Sensory ecology of Ostariophysan alarm substances. Preprints.Org. <https://doi.org/10.20944/preprints201803.0279.v1>

PROTOCOL STATUS











Working

We use this protocol in our group and it is working

GUIDELINES

Guidelines for humane euthanasia must be followed, and animals should NEVER be euthanized without deep anesthesia. See Matthews & Varga (doi: [10.1093/ilar.53.2.192](https://doi.org/10.1093/ilar.53.2.192)) for a review.

MATERIALS

NAME	CATALOG #	
 Distilled Water		Contributed by users
 15 ml sterile falcon tubes and rack		Contributed by users
 scalpel blades		Contributed by users
 Tweezers		Contributed by users
 Ice		Contributed by users
 Glass Petri dishes 90 x 15 cm		Contributed by users
 Graduated Plastic Beakers, 2000ml	BK044-2K.SIZE1	by Bio Basic Inc.
 Cotton Gauze Swabs (10 x 10cm x 8 ply)	G5201 Emh40/0404	Contributed by users
 Disposable pasteur pipettes	EA61.1	by Carl Roth
 Graduated measuring cylinder 10 ml	View	Contributed by users

NAME	CATALOG #
 Aquarium thermometer	View

Contributed by users

SAFETY WARNINGS

1. Alarm substance extraction begins by cutting the fish skin with a scalpel; therefore, care must be taken to avoid accidents
2. All manipulations must be made using surgical gloves, to reduce the risk of contamination by aquatic bacteria (e.g., *Mycobacterium marinum*)

Preparation

- 1 Separate a donor fish from the home tank and transport to the experiment area

Anesthesia

- 2 Make crushed ice using cubes made from mineral water



REAGENT

Ice

Contributed by users

- 3 Fill a 2L beaker with water from the fish facility



REAGENT

Graduated Plastic
Beakers, 2000ml

by [Bio Basic Inc.](#)

Catalog #: [BK044-2K.SIZE.1](#)

- 4 Bring the water temperature to 17 °C by adding the ice chips

- 5 Use a net to transfer the fish from the container to the cold water

- 6 Slowly add ice chips to the beaker to bring the temperature down to 12 °C. The animal is considered anesthetized when it no longer reacts to gentle manipulation.

Alarm substance preparation

- 7 Transfer the anesthetized fish to a Petri dish



REAGENT

Glass Petri dishes 90 x 15

cm

Contributed by users

8 Quickly decapitate the animal using a scalpel; discard the head

9 Remove excess blood using a swab



REAGENT

Cotton Gauze Swabs (10 x

10cm x 8 ply)

Contributed by users

Catalog #: [G5201 Emh40/0404](#)

10 Transfer the body to another Petri dish

11 Using a scalpel, make 15 superficial cuts in both sides of the trunk



REAGENT

scalpel blades

Contributed by users

12 Wash the body with 10 ml ddH₂O



REAGENT

Distilled Water

Contributed by users

13 Discard the body

14 Remove scales and other impurities from the solution

15 Transfer 7 ml from the solution to a Falcon tube

**REAGENT**

15 ml sterile falcon tubes
and rack

Contributed by users

16 Store on ice or refrigerator for up to 4 h.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited