

Preparation PCRS11-Red Sea Medium

The RCC Team

Abstract

The RCC Team, Station Biologique de Roscoff, CNRS/Université P. & M. Curie, Place G. Teissier, 29680 Roscoff, France

Citation: The RCC Team Preparation PCRS11-Red Sea Medium. **protocols.io**

dx.doi.org/10.17504/protocols.io.nusdewe

Published: 15 Mar 2018

Protocol

Preparation of Hepes-NaOH 1M Stock Solution (skip if already available)

Step 1.

To 250mL of H₂O, add gradually 119.15g of Hepes. Adjust pH at 7.5 and complete the volume at 500mL. Store in refrigerator.

 **AMOUNT**

500 ml Additional info: Water

 **AMOUNT**

119.15 g Additional info: Hepes

Preparation of Na₂-EDTA/FeCl₃ Stock Solution (skip if already available)

Step 2.

- To 40mL of HCl 0.1N, add gradually 1.080g of FeCl₃.

- To 40mL of NaOH 0.1N, add gradually 1.488g of Na₂-EDTA.

Now mix both solutions and fill up to the final volume of 2L with sterile water. Store in refrigerator

 **AMOUNT**

1.08 g Additional info: FeCl₃

 **AMOUNT**

40 ml Additional info: HCl 0.1N

 **AMOUNT**

40 ml Additional info: NaOH 0.1N

 **AMOUNT**

1.448 g Additional info: Na₂-EDTA

Preparation of Sodium Phosphate Stock Solution (skip if already available)

Step 3.

Prepare two solutions:

- Monosodium dihydrogen phosphate (NaH_2PO_4) at 50mM (6g in 1L)
- Disodium hydrogen phosphate (Na_2HPO_4) at 50mM (3.55g in 500mL)

Make an equimolar mixture of this two solutions and adjust the pH at 7.5.

☐ CONCENTRATION

0.05 Molarity (M) Additional info: NaH_2PO_4

☐ CONCENTRATION

0.05 Molarity (M) Additional info: Na_2HPO_4

Preparation of Trace metals "Gaffron+Se" Stock Solution (skip if already available)

Step 4.

To 500mL of H_2O , add gradually the following nutrients

Quantity (mg/L)	Compound Final	Concentration in media (nM)
186	Boric acid (H_3BO_3)	150
101	Manganese (II) Sulfate Monohydrate ($\text{MnSO}_4 \cdot \text{H}_2\text{O}$)	30
1.98	Sodium Tungstate dihydrate ($\text{Na}_2\text{WO}_4 \cdot 2\text{H}_2\text{O}$)	0.3
5.16	Ammonium molybdate tetrahydrate ($(\text{NH}_4)_6\text{Mo}_7\text{O}_{24} \cdot 4\text{H}_2\text{O}$)	1.45
7.14	Potassium bromide (KBr)	3
4.98	Potassium iodide (KI)	1.5
17.25	Zinc sulfate heptahydrate ($\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$)	3
9.25	Cadium Nitrate ($\text{Cd}(\text{NO}_3)_2 \cdot 4\text{H}_2\text{O}$)	1.5
8.76	Cobalt (II) Nitrate ($\text{Co}(\text{NO}_3)_2 \cdot 6\text{H}_2\text{O}$)	1.5
7.5	Copper (II) Sulfate ($\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$)	1.5
7.1	Nickel Chloride ($\text{NiCl}_2 \cdot 6\text{H}_2\text{O}$)	1.5
2.4	Chromium (III) Nitrate ($\text{Cr}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$)	0.3
1.5	Vanadyl Sulfate Pentahydrate ($\text{VOSO}_4 \cdot 5\text{H}_2\text{O}$)	0.3
28.4	Aluminium Potassium Sulfate ($\text{KAl}(\text{SO}_4)_2 \cdot 12\text{H}_2\text{O}$)	3
3.3	Selenium (IV) Oxyde (SeO_2)	1.5

Complete the volume at 1L and store in refrigerator.

Preparation PCRS11-Red Sea medium

Step 5.

To 1L of H_2O , add 33,33g of Red Sea Salt. Dissolve by shake (20min on agitator)

Preparation PCRS11-Red Sea medium

Step 6.

Heat seawater during 20min at 100°C.

TEMPERATURE

100 °C Additional info:

Preparation PCRS11-Red Sea medium

Step 7.

Under the hood, to water, add these nutrients beforehand autoclaved (except vitamin):

Quantity	Compound	Final Concentration
1.0 mL	Hepes-NaOH 1M (pH 7.5) (See recipe below)	1mM
1.0 mL	Na ₂ -EDTA/FeCl ₃ (See recipe above)	8μM
1.0 mL	Sodium Phosphate (NaPO ₄) 50mM (pH 7,5) (See recipe above)	50 μM
1.0 mL	Ammonium Sulfate 400mM (NH ₄) ₂ -SO ₄	400μM
0,1 mL	Trace metals “Gaffron+Se” (See recipe above)	
0.1 mL	Cyanocobalamin 10mg/L (Vit. B12)	1μg/L

Filter the medium on 0.2 microns.