

Jun 18,  
2019

## Phagocytosis Bead Conjugation V.2

Andrew Crowley<sup>1</sup>

<sup>1</sup>Dartmouth College, Ackerman Lab

Working

dx.doi.org/10.17504/protocols.io.4cygsxw

Ackerman Lab



Andrew Crowley

Dartmouth College, Ackerman Lab



### ABSTRACT

Protocols for the amine coupling of fluorescent microspheres with antigen or capture reagent for opsonization and phagocytosis

Quantity as written: approx.  $1 \times 10^8$  under typical yields (enough for 400 phagocytosis wells under standard protocol)

### STEPS MATERIALS

NAME	CATALOG #	VENDOR
PBS-TBN	P0210	Teknova
PBS-TBN	P0210	Teknova

### Activation

- 1 Calculate the beads / mL from the percent solids as follows:

$$(6W \times 10^{12}) / (\rho \times \pi \times \phi^3)$$

W = grams solid per mL

$\rho$  = density (g/mL)

$\phi$  = particle diameter ( $\mu\text{m}$ )

polystyrene density = 1.055 /mL

- 2 Resuspend beads by briefly vortexing; dispense **200  $\mu\text{l}$**  (contains approx. 7 billion beads)

- 3 Pellet by centrifugation **17000 x g** **00:03:00** and remove supernatant by pipette

- 4 Resuspend in a total of **1 mL** consisting of:

- sulfo-NHS **50 mg/mL** made in **10 Millimolar (mM)** MES (pH 5.0) **100  $\mu\text{l}$**

- EDC **50 mg/mL**, made in **10 Millimolar (mM)** MES (pH 5.0) **100  $\mu\text{l}$**

- **10 Millimolar (mM)** MES (pH 5.0) **800  $\mu\text{l}$**

5 Incubate with end-over-end mixing (shielded from light) ⌚ 00:20:00 🌡 Room temperature

#### Coupling

6 Wash by pelleting beads and replacing supernatant with [M]10 Milimolar (mM) MES (pH 5.0) 🧴 1 ml  
🌀 17000 x g ⌚ 00:02:00

7 Formulate protein for coupling @ approx. [M]25 µg/mL in [M]10 Milimolar (mM) sodium acetate (pH 5.0) 🧴 7.5 ml



n.b.

The large volume used for the step is based on observations that sufficiently high bead concentrations during this step may lead to crosslinking between microspheres (**anecdotal**)

8 Incubate with end-over-end mixing (shielded from light) ⌚ 03:00:00 🌡 Room temperature

#### Blocking + Washes

9 Collect the full volume in the base of the reaction tube by centrifugation 🌀 250 x g ⌚ 00:01:00

10 Split the 🧴 7.5 ml volume into multiple fractions in microcentrifuge tubes (🧴 2 ml) in order to centrifuge at a sufficiently high speed  
🌀 17000 x g ⌚ 00:03:00

11 Remove supernatant from each fraction by pipette

12 Resuspend **one** fraction in PBSF (ie. 1x PBS + [M]0.1 Mass/Volume Percent BSA)  1 ml

13 Transfer the entire contents of the first fraction to the second and resuspend; continue through the fractions until as many beads as possible have been collected in a single microcentrifuge tube

14 Pellet beads via centrifuge  17000 x g  00:02:00

15  go to step #11 for a total of 3 washes

16 Resuspend in PBSF  500 µl


17 Will the beads be used for an assay on the same day as production?

\_\_\_\_\_ step case \_\_\_\_\_

No

no description provided



18 Block via  16:00:00 (ie. overnight)  4 °C

#### Counting + Qualification

19 Make serial dilutions of 1:100, 1:1000, and 1:10,000

20 Measure concentration by flow cytometry using the fluorescent dye of the microsphere as the trigger condition




Typical yield is approx.  $1 - 2 \times 10^8$  beads, or around 25% of the original amount



MACS Quant  
flow cytometer  
Miltenyi 130-096-343

21 In a non-binding, 96-well plate, formulate serial dilutions of positive and negative control antibodies beginning at [M]5 µg/mL and diluted by 1:3 over an approx. 4 point series

22 Add 100,000 beads to each well

- 23 Incubate with shaking (shielded from light) ⌚ 01:00:00 🌡 Room temperature
- 24 Pellet by centrifugation 🌀 3500 x g ⌚ 00:05:00 and remove supernatant by decanting
- 25 Resuspend beads in PBST (ie. 1x PBS + [M] 0.05 Volume Percent Tween20) 📄 200 µl
- 26 [go to step #24](#) for a total of 2 washes
- 27 Resuspend beads in fluorescent secondary 📄 100 µl [M] 0.65 µg/mL
-  n.b.

Make sure that the fluorescent molecule chosen for the secondary does not have spectral overlap with the dye in the microspheres themselves
- 28 Incubate with shaking (shielded from light) ⌚ 00:30:00 🌡 Room temperature
- 29 Pellet by centrifugation 🌀 3500 x g ⌚ 00:05:00 and remove supernatant by decanting
- 30 Resuspend beads in PBST (ie. 1x PBS + [M] 0.05 Volume Percent Tween20) 📄 200 µl
- 31 Pellet by centrifugation 🌀 3500 x g ⌚ 00:05:00 and remove supernatant by decanting
- 32 Resuspend beads in 1x PBS 📄 100 µl
- 33 Measure fluorescent intensity by flow cytometry



MACS Quant  
flow cytometer  
Miltenyi 130-096-343

#### Washing + Blocking

step case

No

no description provided

NaN Block via 🕒 16:00:00 (ie. overnight) 🌡 4 °C

#### Counting + Qualification

NaN

#### Washing + Blocking

step case

Yes

no description provided

NaN Block via 🕒 00:30:00 🌡 Room temperature

#### Blocking + Washes

step case


Yes

no description provided

18 Block via 🕒 00:30:00 🌡 Room temperature

#### Counting + Qualification

19 Make serial dilutions of 1:100, 1:1000, and 1:10,000



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited