



Apr 07, 2019

In devel.

## RNA isolation for tissue

Version 2

Chin Yee Tan<sup>1</sup><sup>1</sup>Duke University[dx.doi.org/10.17504/protocols.io.zvdf626](https://doi.org/10.17504/protocols.io.zvdf626)

- 1 Homogenization in TriZol  
Add 1 ml of Trizol reagent per 30 mg of tissue and homogenize using handheld homogenizer
- 2 Incubate at RT for 00:05:00 to allow nucleoprotein complexes to dissociate
- 3 Add 1/5 the volume of Trizol ( 0.2 ml ) of chloroform carefully, and vortex to mix well
- 4 Spin down at max speed in a chilled centrifuge for 00:15:00
- 5 Carefully remove the top aqueous phase and transfer to a new Eppendorf tube  
The interphase and bottom organic phase can be saved for DNA and protein respectively
- 6 To the aqueous phase, add 500 µl of 100% isopropanol, mix by inversion and incubate at -20 °C for a minimum for 02:00:00
- 7 Spin down at max speed for 00:30:00 s to precipitate RNA
- 8 Remove supernatant, and add 1 ml of 75% Ethanol to wash the pellet
- 9 Spin down at max speed for 00:15:00 and remove supernatant
- 10 Resuspend pellet in appropriate volume of nuclease free water

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits

