



Aug 20,
2019

RNA Isolation from Plant Tissue Protocol 5: pBIOZOL Method

1

Works for me

dx.doi.org/10.17504/protocols.io.4rdgv26



Scott C. Edmunds
GigaScience/BGI Hong Kong/Bauhinia Genome



ABSTRACT

Implemented by: Beijing Genomics Institute

This protocol is part of a collection of eighteen protocols used to isolate total RNA from plant tissue. (RNA Isolation from Plant Tissue Collection: <https://www.protocols.io/view/rna-isolation-from-plant-tissue-439gyr6>)

journal.pone.0050226.s0
11.PDF

MATERIALS TEXT

Reagents

- 5 M NaCl
- Chloroform
- Isopropyl alcohol
- 75 % ethanol (DEPC treated)
- pBIOZOL Reagent (Beijing Bai billion New Technology Co., Beijing, China)
- RNase-free water

SAFETY WARNINGS

Please see SDS (Safety Data Sheet) for hazards and safety warnings.

- 1 Grind tissue to a powder in liquid nitrogen.
- 2 Add 1.3 ml of cold (4 °C) pBIOZOL Reagent for up to 100 mg of frozen, ground tissue.
 - 2.1 Mix by briefly vortexing or flicking the bottom of the tube until the sample is thoroughly suspended.
- 3 Incubate the tube for 00:05:00 at 20 °C.



Lay the tube down horizontally to maximize surface area during RNA extraction.

- 4 Centrifuge for 🕒 00:02:00 at 🌀 12000 x g.
- 4.1 Transfer the supernatant to an RNase-free tube.
- 5 Add 📏 100 µl of [M]5 Molarity (M) NaCl to the extract.
- 5.1 Tap tube to mix.
- 6 Add 📏 300 µl of chloroform.
- 6.1 Mix thoroughly by inversion.
- 7 Centrifuge the sample at 🌡 4 °C for 🕒 00:10:00 at 🌀 12000 x g to separate the phases.
- 7.1 Transfer the top aqueous phase to an RNase-free tube.
- 8 Add to the aqueous phase an equal volume of isopropyl alcohol.
- 8.1 Mix.
- 8.2 Let stand at 🌡 20 °C for 🕒 00:10:00 .
- 9 Centrifuge the sample at 🌡 4 °C for 🕒 00:10:00 at 🌀 12000 x g.
- 10 Decant the supernatant, taking care not to lose the pellet.

10.1 Add  1 ml of chilled 75 % ethanol to the pellet.



Pellet may be difficult to see.

11 Centrifuge at  Room temperature for  00:05:00 at  12000 x g.

11.1 Decant the liquid carefully, taking care not to lose the pellet.

11.2 Briefly centrifuge to collect the residual liquid and remove it with a pipette.

12 Add  10 µl –  30 µl RNase-free water to dissolve the RNA.

12.1 Pipette the water up and down over the pellet to dissolve the RNA.



If any cloudiness is observed, centrifuge the solution at  Room temperature for  00:01:00 at  12000 x g and transfer the supernatant to a fresh tube.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited