

Apr 24, 2019

Working


Dot blot to verify microinjection efficacy in *Euplotes crassus*

Forked from [Dot blot to verify microinjection efficacy in *Euplotes crassus*](#)RACHELE CESARONI¹, [Rachele Cesaroni](#)²¹University of Bern, Institute of Biology, ²Universität Bern[dx.doi.org/10.17504/protocols.io.2aqqadw](https://doi.org/10.17504/protocols.io.2aqqadw)

Protist Research to Optimize Tools in Genetics (PROT-G)

Angela Piersanti
University of Camerino

- 1 Grow cells up to a density of 1000 cells/ml, and transfer 400 µl of each culture to an Eppendorf tube.
- 2 Pellet the cells at max speed for 5 minutes, and resuspend them in 400 µl of ddH₂O.
- 3 Add 50 µl of 0.5 M EDTA, pH 8.0 and 50 µl of 4 M NaOH to the cells.
- 4 Lyse the cells at 68°C for 30 minutes, and centrifuge them for 1 minute at 13000 rpm, then chill them on ice.
- 5 Soak a membrane for nucleic acid binding and a sheet of Whatmann paper in 0.4 M of NaOH, and place both in a Dot blot apparatus.
- 6 Switch on the vacuum pump, and transfer all the lysates into the chamber.
- 7 Soak the membrane in 0.4 M NaOH for 15 minutes, and wash it in 2X SSC buffer for 10 minutes, then drain it on a sheet of Whatmann paper.
- 8 Pre-hybridize the membrane in 20 ml of Church buffer at 60°C for 2 hours, then hybridize it at constant temperature (annealing temperature) o/n in 10 ml of Church buffer containing the radioactive probe.
- 9 Wash the membrane twice with 2X SSC, 0.1% SDS for 15 minutes each time.
- 10 Place the membrane in a plastic foil, expose it overnight and scan it to determine which clones have sufficient DNA from the microinjection (compared to suitable DNA standards used as positive controls).



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited