

100x BG-11 media

Anna Behle, Miriam Dreesbach, Susanne Vollmer

Abstract

This protocol describes how to prepare a 100x BG-11 +N stock.

1x BG-11 media could be created from the stock and supplemental stocks and trace metals could be added afterwards.

This protocol is based on Anne Behle M.Sc. [Recipe for standard BG-11 media](#) protocol.

Recipes for standard and alternative BG11 for culturing freshwater cyanobacteria, such as *Synechocystis* sp. PCC 6803, as described.

Media is usually not suitable for marine cyanobacteria.

Citation: Anna Behle, Miriam Dreesbach, Susanne Vollmer 100x BG-11 media. **protocols.io**

dx.doi.org/10.17504/protocols.io.n97dh9n

Published: 06 Apr 2018

Guidelines

Always work under sterile conditions when handling sterile media or stocks. Work under the clean bench.

Protocol

Safety first

Step 1.

Always work under sterile conditions.

Step 2.

Add $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ ($3.6 \text{ g} \cdot \text{L}^{-1}$). We did not have hydrated CaCl_2 . Therefore, we used 2.718 g CaCl_2 (powder) for 1L of the stock.

Step 3.

Add Citric acid ($0.6 \text{ g} \cdot \text{L}^{-1}$)

Step 4.

Add NaNO_3 ($149.58 \text{ g} \cdot \text{L}^{-1}$)

Step 5.

Add $\text{MgSO}_4 \cdot 7 \text{ H}_2\text{O}$ ($7.49 \text{ g} \cdot \text{L}^{-1}$). We do not have hydrated MgSO_4 . Therefore, we used 3.659 g MgSO_4 (powder) for 1L of the stock.

Step 6.

Add 0.25 M $\text{Na}_2\text{-EDTA}$, pH 8.0 ($0.56 \text{ ml} \cdot \text{L}^{-1}$). We did not have a 0.25 M $\text{Na}_2\text{-EDTA}$ stock. Therefore, we dissolved 2.32 g $\text{Na}_2\text{-EDTA}$ (powder) and filled the bottle up with MilliQ (nuclease-free) water up to 25 mL . Afterwards, the pH was adjusted with diluted NaOH to $\text{pH} = 8.0$.

Step 7.

The bottle was filled with MilliQ (nuclease-free) water up to 1 L and mixed a big stir bar until the powder was completely dissolved and the solution was clear.

Step 8.

Afterwards, the stir bar was removed and the media was directly autoclaved at $120 \text{ }^\circ\text{C}$ for 30 Minutes.

Warnings

Wear gloves when preparing stocks!