

# Viral metagenomic analysis on Cabbage Patch Kids

Stanton Burnton

## Abstract

This protocol is a parody of an existing scientific journal article and is intended to showcase the features available on the protocols.io platform in order for VERVE Net members to develop future protocols. While the viral metagenomic analysis on Cabbage Patch Kids is impossible to perform, due to the real-world non-existent nature of Cabbage Patch kids, the sourced article is *Metagenomic analysis of viruses associated with field-grown and retail lettuce identifies human and animal viruses* and the "Material and methods" section has been adapted from using actual lettuce data to a metagenomic protocol for Cabbage Patch Kids.

**Citation:** Stanton Burnton Viral metagenomic analysis on Cabbage Patch Kids. **protocols.io**

dx.doi.org/10.17504/protocols.io.eqtbdwn

**Published:** 28 Mar 2016

## Before start

It is advisable to obtain all necessary materials prior to obtaining field Cabbage Patch Kids. Cabbage Patch Kids mutate quickly from field removal and must be maintained in an adequate environment to prevent a syndrome what has recently been described as *Ravage Patch Kids*.

## Protocol

### Cabbage Patch Kid (CPK) collection for metagenomic analysis

#### Step 1.

Collect Cabbage Patch Kids (CPK) at varying stages of farming production, including:

- hand-cut by field workers and packaged
- washed and cut by processing workers
- hand-cut by researchers using sterile equipment
- produce distribution centers

### Cabbage Patch Kid (CPK) collection for metagenomic analysis

#### Step 2.

Use sterile gloves and cut outer leaflets off using a scalpel, careful not to harm the body of the CPK, before placing CPK in large sterile Whirl-pak bags.

### Virus recovery from CPK

#### Step 3.

Wash each sample in the Whirl-pak bag with 250 ml sterile 100 mM Tris — 50 mM glycine buffer at a pH of 9.5 and gently mix for 20 min at room temperature. Recover the wash solution immediately and adjust the pH to neutral  $7.2 \pm 0.2$

📄 **AMOUNT**

250 ml Additional info:

## DURATION

00:20:00

## PROTOCOL

### . [Tris-glycine buffer](#)

CONTACT: [Stanton Burnton](#)

#### Step 3.1.

Mix 100mM Tris with 50mM glycine to yield a buffer

#### AMOUNT

250 ml Additional info:

#### Virus recovery from CPK

#### Step 4.

Use polyethylene glycol (PEG) precipitation to concentrate and purify the viral particles contained in the wash solution.

#### Virus recovery from CPK

#### Step 5.

Mix samples with 10% (weight/volume) PEG 8000 and .3 M NaCl (weight/volume), then incubate at 4 °C for 18 hours before centrifuging the samples at 10,800 × *g* (8000 rpm) for 30 min at 4 °C.

#### DURATION

18:30:00

#### Virus recovery from CPK

#### Step 6.

Pour off the supernatant and dissolve the pellet in 20 mL of sterile phosphate buffered saline, letting soak for 1 hour at room temperature

#### AMOUNT

20 ml Additional info:

#### REAGENTS

✓ 1X PBS (Phosphate-buffered saline ) by Contributed by users

#### DURATION

01:00:00

#### Virus recovery from CPK

#### Step 7.

Add an equal volume of chloroform to each PEG precipitate to remove the PEG and purify the sample.

#### REAGENTS

✓ Chloroform by Contributed by users

#### Virus recovery from CPK

#### Step 8.

Vortex the solutions for 30 seconds and centrifuge at 3000 × *g* (4300 rpm) for 15 min at 4 °C to collect the supernatant containing virus particles.

#### DURATION

00:15:30

#### Virus recovery from CPK

#### Step 9.

Pass the remaining supernatant through 0.45 and 0.22 µm filters and further concentrate to approximately 1 mL by Amicon centrifugal ultrafiltration (30 kDa)

## Nucleic acid extraction and sequencing

### Step 10.

Treat the final 1 mL concentrates with 100 units of DNase-I for 1 hour at 37 °C before nucleic acid extraction to remove free nucleic acids from the concentrated virus samples.



#### REAGENTS

deoxyribonuclease I (DNase I, 100 U/ml) D 4263 by [Sigma Aldrich](#)



#### DURATION

01:00:00

## Nucleic acid extraction and sequencing

### Step 11.

Extract viral DNA and RNA using a PureLink viral RNA/DNA mini kit (Life Technologies) following the manufacturer's instructions. For each viral concentrate, prepare three individual nucleic acid extracts to minimize nucleic acid extraction bias.

## Nucleic acid extraction and sequencing

### Step 12.

Following extraction, screen the samples with 16S ribosomal DNA (rDNA) PCR with 27F/1492R universal primers to ensure the absence of any contaminating microbial DNA.



#### NOTES

**Stanton Burnton** 21 Mar 2016

To obtain a sufficient DNA and cDNA (for RNA viruses) for metagenomics sequencing, the viral nucleic acids were reverse transcribed and amplified as previously described

## Nucleic acid extraction and sequencing

### Step 13.

Reverse transcribe RNA with Primer A (5'-GTTTCCCAGTCACGATCNNNNNNNNN-3') using Superscript III reverse transcriptase

## Nucleic acid extraction and sequencing

### Step 14.

Use Sequenase 2.0 for second-strand cDNA synthesis and for random-primed amplification of viral DNA. Subject each sample to 40 cycles of PCR amplification with Primer B (5'-GTTTCCCAGTCACGATC-3') using AmpliTaq Gold

## Nucleic acid extraction and sequencing

### Step 15.

Perform three PCR reactions from the same nucleic acid extract to minimize amplification bias and pool the PCR products. Purify PCR products using Promega Wizard SV Gel and a PCR Clean-Up System

## Nucleic acid extraction and sequencing

### Step 16.

Prepare libraries from each sample using a Rubicon ThruPLEX DNA-seq kit with a unique dual index adapter pair for each sample. Sequence samples in a 2 × 100-base pair (bp) paired end format using two lanes of an Illumina HiSeq 2500 Rapid Run flow cell

## Bioinformatics analysis of viromes

### Step 17.

Screen each dataset for the 17-bp Primer B sequence and any reads homologous to the Primer B sequence at their 5' ends, removing using cutadapt with a maximum error rate of 0.2 and minimum overlap of 10 bases



#### NOTES

**Stanton Burnton** 21 Mar 2016

Prior to assembly and annotation of the metagenomic dataset, the quality of the Illumina sequencing data was checked using FastQC

#### Bioinformatics analysis of viromes

##### Step 18.

Use Trimmomatic for sequencing adapter removal and quality trimming with parameters including : a maximum mismatch count value of 2 allowed for a full match (seed mismatch), a palindrome clip threshold of 30, a simple clip threshold of 10, a minimum adapter length of 8 with both the forward and reverse read kept, removal of low quality leading and trailing bases below a quality of 3, a 4-base sliding window scan that cuts when the average quality is below 15, and removal of reads less than 30 bases long

#### Bioinformatics analysis of viromes

##### Step 19.

Following filtering and trimming of raw reads, subject paired-end reads to *de novo* assembly into a longer contiguous sequence (contig) using IDBA-UD

#### Bioinformatics analysis of viromes

##### Step 20.

Query Contigs larger than 200 bp were then against the National Center for Biotechnology Information (NCBI) Viral Reference Sequence (RefSeq) database for taxonomic assignment using BLASTX with an E-value cutoff of  $10^{-5}$

#### Bioinformatics analysis of viromes

##### Step 21.

Parse the BLASTX output using the MEtaGenome Analyzer (MEGAN) version 5.6.6 with the following parameters for the Lowest Common Ancestor (LCA) algorithm: min score = 50.0, max expected =  $1.0 \times 10^{-5}$ , top percent = 10.0, min support percent = 0.1, min support = 1, and LCA percent = 100.0

#### Bioinformatics analysis of viromes

##### Step 22.

Extract contigs identified as viral pathogens of human and animal and use them as the queries in BLASTX against the NCBI non-redundant (nr) sequence database.

#### 🔗 NOTES

**Stanton Burnton** 21 Mar 2016

For rotavirus contigs, BLASTN was used

#### Bioinformatics analysis of viromes

##### Step 23.

To determine relative abundance of a phylogenetic group, perform read mapping to contigs using Bowtie 2 version 2.1.0 with default settings

#### Bioinformatics analysis of viromes

##### Step 24.

Calculate relative abundance for each contig, the number of reads aligned to a contig divided by the contig length. Calculate the relative abundance of each phylogenetic group by summing the abundance of each contig classified in a particular group.

## Warnings

Take care not to harvest live Cabbage Patch Kids, as legislation prohibits the possession, distribution, and sale of live Cabbage Patch Kids. Please refer to the government guidelines set in place which detail the legalities of Cabbage Patch Kid farming and production.

**This protocol is derived from a journal article in the International Journal of Food Microbiology about lettuce and not actually a protocol used to conduct a viral metagenomic analysis on Cabbage Patch Kids.**