

Making Differentiation Media for SH-SY5Y

Kenneth Schackart¹

¹University of Arizona

[dx.doi.org/10.17504/protocols.io.xuufnww](https://doi.org/10.17504/protocols.io.xuufnww)

 Kenneth Schackart 

ABSTRACT

This describes how to make differentiation medium for SH-SY5Y using all-*trans*-retinoic acid.

PROTOCOL STATUS




In development




We are still developing and optimizing this protocol

MATERIALS TEXT

- all-*trans*-retinoic acid
- B-27 Supplement
- L-glutamine
- Neurobasal Medium (NBM)
- 100% Ethanol
- Autoclaved deionized water (ADIW)

Make ATRA Stock Solution

- 1 In a 15 mL centrifuge tube, add  14.25 ml ethanol and  0.75 ml autoclaved deionized water. This will make a 95% ethanol solution.
- 2 Transfer  3330 µl of 95% EtOH solution to another 15 mL centrifuge tube.
- 3 Weigh out  5 mg ATRA.

 **NOTE**
ATRA is light sensitive, so do this step with the lights off.
- 4 Add the  5 mg ATRA to the  3330 µl 95% EtOH solution. This makes a 5 mM stock ATRA solution.
- 5 Vortex mix the solution.

- 6 This solution is used at a rate of stock solution per differentiation media. The final concentration should be 10 μM.

NOTE

We most recently made 10 aliquots, each with . This is to make 10 media batches, each using .

- 7 Stock solution can be stored for up to 6 weeks at .

Make l-Glutamine Stock Solution

- 8 Weigh out l-glutamine.
- 9 Aliquot NBM into a 15 mL centrifuge tube.
- 10 Add the l-glutamine to the NBM.
- 11 Vortex mix the solution.
- 12 Aliquot the solution into aliquots.
This solution is 20x. The concentration is 80 mM and the final desired concentration is 4 mM.

NOTE

This protocol makes 10 aliquots of . Each stock aliquot makes differentiation medium.

- 13 Store aliquots at .

Aliquot B-27 Supplement

- 14 Aliquot batches of 50x B-27 supplement.
- 15 Each aliquot will make differentiation media.

16 Store aliquots at  -20 °C .

Making the media

17





NOTE

Do not mix the media until you are ready to change the media in a flask. It degrades quickly.

Warm  30 ml NBM in water bath.

18 Thaw one 50x B-27 supplement aliquot, and one 20x l-glutamine supplement.

19 In biosafety cabinet, combine the following in a 50 mL tube

-  29.34 ml Nuerobasal Medium
-  600 µl 50x B-27 Supplement
-  1.5 ml 20x l-glutamine Stock Solution
-  60 µl 500x All-*trans*-retinoic Acid Stock Solution

20 Sterile filter the differentiation media.

21 Use immediately.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited