

Oct 18, 2019

Frequently sampled Insulin glucose tolerance test V.2 👄

Timothy Nichols¹, David Clemmons¹

¹University of North Carolina at Chapel Hill



Diabetic Complications Consortium Tech. support email: rmcindoe@augusta.edu



ABSTRACT

Summary:

This assay is used by the DiaComp to measure glucose tolerance and insulin sensitivity in pigs.

Diabetic Complications:



Reference:

1. Bergman RN, Finegood DT, Ader M: Assessment of insulin sensitivity in vivo. Endocr Rev 1985, 6:45-86.

EXTERNAL LINK

https://www.diacomp.org/shared/document.aspx?id=50&docType=Protocol

MATERIALS

NAME CATALOG # VENDOR YSI Glucose Analyzer **YSI Life Sciences**

ICN Insulin RIA kit

MATERIALS TEXT

Reagents Quantity Required

Reagent/Material	Quantity Required
Intravenous catheter	2
ICN Insulin RIA kit	1 kit
YSI Glucose Analyzer	1

1 FSIGT or Bergman analysis

Pigs are studied after an overnight fast. The food intake of the animals is monitored for 3 days prior to the fast to ensure adequate carbohydrate intake. Two intravenous catheters are placed, one for sampling and one for infusing glucose and insulin. A bolus of glucose (0.3 gm/kg) is administered as a 50% solution over \sim 5 min. Blood samples are obtained at \sim 15, \sim 10, \sim 5, \sim 1, 0, 2, 3, 4, 5, 6, 8, 10, 12, 14, 16, and 19 minutes. At 20 minutes an insulin bolus (0.03U/kg) is injected and frequent blood samples for insulin and glucose measurements are collected up to the 180 minute time point. Insulin is measured by RIA (ICN) and glucose is measured on a YSI instrument (Yellow Springs, Ohio). The data were analyzed by the Bergman method to calculate an insulin sensitivity index (S_I) using MINMOD Millennium version 6.02.1

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited