



Prepare 2L of B-Broth for culturing bacteria

Cancer Research UK / Wellcome Gurdon Institute media kitchen¹

¹Cancer Research UK / Wellcome Gurdon Institute

[dx.doi.org/10.17504/protocols.io.3nmgmc6](https://doi.org/10.17504/protocols.io.3nmgmc6)

Working



Cristian Riccio ⚡

ABSTRACT

Prepare 2L of B-Broth for culturing bacteria

MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
double distilled water (ddH2O)		
Sodium chloride meets analytical specification of Ph.Eur Fisher Chemical	S/3160/65	Fisher Scientific
bacto tryptone	View	

SAFETY WARNINGS

Make sure your know how to use the autoclave.

1 B-Broth.xls

2

Ingredients	Quantity	
Bacto tryptone	20g	
NaCl		10g
Double distilled H2O	up to 2L	

3

1	Dissolve Bacto tryptone and NaCl in 1L double distilled H2O in a 2L measuring cylinder
2	Make up to final volume using double distilled H2O
3	Aliquot in 100 ml bottles, label and autoclave



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited