

Version 3

Oct 17, 2018

Working

Step 2: Preparing amino acid, polyphosphates, and maltodextrin-based energy solutions for cell-free reactions Version 3

Anibal Arce Medina¹, Juan Puig Contreras², Fernan Federici²¹Pontificia Universidad Católica de Chile, ²Pontificia Universidad Católica de Chile

dx.doi.org/10.17504/protocols.io.urzev76

SynBioUC Chile



Anibal Arce Medina

Pontificia Universidad Catolic...



















PROTOCOL STATUS

Working

We use this protocol in our group and it is working

STEPS MATERIALS

NAME	CATALOG #	VENDOR
 L-aminoacids	LAA21-1KT	Sigma Aldrich
 L-Proline	81709-25G	Sigma Aldrich
 L-Cysteine	30089	Sigma Aldrich
 L-histidine	53319-25G	Contributed by users
 L-lysine	L5501-5G	Contributed by users
 DTT	D0632	Sigma Aldrich
 CoA	C4282	Sigma Aldrich
 Folinic Acid	F7878	Sigma Aldrich
 Spermidine	85558	Sigma Aldrich
 NAD (β-Nicotinamide adenine dinucleotide hydrate)	N6522	Sigma Aldrich
 HEPES	H6147	Sigma Aldrich
 ATP (Adenosine 5'-triphosphate dipotassium salt hydrate)	A8937	Sigma Aldrich
 GTP (Guanosine 5'-triphosphate sodium salt hydrate)	10106399001 ROCHE	Sigma Aldrich
 UTP (Uridine 5'-triphosphate trisodium salt dihydrate)	94370	Sigma Aldrich
 Maltodextrin (4.0-7.0 dextrose equivalent)	419672	Sigma Aldrich
 sodium hexametaphosphate	305553	Sigma Aldrich

Prepare individual amino acid stock solutions

1

Dissolve the given amount of each amino acid with 1 ml of 5M KOH in a 3 ml screw-cap tube according to the following list.

Aminoacid	Amount of powder needed
Alanine	445.45 mg

Arginine	644.6502 mg
Asparagine	660.6 mg
Aspartic acid	665.5 mg
Cysteine	370.7496 mg
Glutamic acid	735.65 mg
Glutamine	447.219 mg
Glycine	375.35 mg
Histidine	775.75 mg
Isoleucine	655.9 mg
Leucine	401.4108 mg
Lysine	447.3414 mg
Methionine	456.5826 mg
Phenylalanine	330.38 mg
Proline	575.65 mg
Serine	529.5 mg
Threonine	595.6 mg
Tryptophan	408.46 mg
Tyrosine	554.4414 mg
Valine	358.479 mg

Vortex to dissolve. Each solution can be stored at -20 °C.

 **REAGENT**
L-aminoacids
by [Sigma Aldrich](#)
Catalog #: [LAA21-1KT](#)

 **REAGENT**
L-Proline
by [Sigma Aldrich](#)
Catalog #: [81709-25G](#)

 **REAGENT**
L-Cysteine
by [Sigma Aldrich](#)
Catalog #: [30089](#)

 **REAGENT**
L-histidine
Contributed by users
Catalog #: [53319-25G](#)

 **REAGENT**
L-lysine
Contributed by users
Catalog #: [L5501-5G](#)

Prepare amino acid mix solution (~17 nM each)

- In a 50 mL Falcon tube add the following volumes of the individual amino acid stock solutions

Aminoacid Stock solution in 5M KOH	Volume needed
------------------------------------	---------------

Alanine	136 ul
Arginine	222 ul
Asparagine	136 ul
Aspartic acid	136 ul
Cysteine	222 ul
Glutamic acid	136 ul
Glutamine	222 ul
Glycine	136 ul
Histidine	136 ul
Isoleucine	136 ul
Leucine	222 ul
Lysine	222 ul
Methionine	222 ul
Phenylalanine	340 ul
Proline	136 ul
Serine	136 ul
Threonine	136 ul
Tryptophan	340 ul
Tyrosine	222 ul
Valine	222 ul
TOTAL	3816 µL

Prepare amino acid mix solution (~14 mM each)

- 3 Add 35.084 µL of sterile water and 1.100 µL of Acetic acid (glacial) to the amino acid mix. Vortex well and aliquot in 2.0 ml Eppendorf tubes. Aliquots can be stored at -80°C

Prepare maltodextrin-based energy solution

- 4 Prepare the following stock solutions:

1 M DTT: Put 2.31 g DTT in a 15 ml Falcon tube and fill with water to 15 ml. Sterilize using a 0.22 µm filter. Aliquot in 1.5 ml Eppendorf tubes. Store at -20 °C for later use.

2M HEPES pH 8: Weight 19.1 g HEPES (MW 238.21). Dissolve with 30ml water. Adjust pH to 8.0 with KOH. Fill with water to 40 ml.

5 mg/ml tRNA solution: Put 30 mg of tRNA in a 1.5 ml Eppendorf tube. Fill with water to 600 µl.

CoA stock solution: Put 30 mg of CoA (MW 767.53) in a 1.5 ml Eppendorf tube and fill with water to 600 µl.

38.3 mM NAD solution: Put 34.83 mg of NAD (MW 663.43) in a 1.5 ml Eppendorf tube, add 27 µl of Tris at 2 M (to bring the solution to pH 8.0). Finally, fill with water to 300 µl.

23 mM cAMP solution: Put 42.80 mg of cAMP (MW 329.22) in a 1.5 ml Eppendorf tube, add 73 µl of Tris at 2 M (to bring the solution to pH 8.0). Finally, fill with water to 200 µl.

40 mM Folinic acid solution: To 20 mg of solid folinic acid calcium salt (MW 511.5), add 1.15 ml of water.

Spermidine stock solution: Briefly warm spermidine (two minutes at 37°C) in order to melt it. Then, put 23.55 µl of spermidine (MW 145.25) into a 1.5 Eppendorf tube and fill with water to 150 µl.

240 mg/ml maltodextrin solution: Put 2.4 g of maltodextrin in a 15 ml Falcon tube. Dissolve, and fill with water to 10 ml.

Nucleotide Mix solution: Put 145 mg of ATP dipotassium salt dihydrate (MW 619.4), 133 mg of GTP disodium salt dihydrate (MW 567.14), 79.4 mg of CTP disodium salt dihydrate (MW 563.16), 82.6 mg of UTP trisodium salt dihydrate (MW 586.12) in a 1.5 Eppendorf tube. Add 353 µl of KOH (15 w/v %). Finally, fill with water to 1.5 ml.



REAGENT

DTT

by [Sigma Aldrich](#)

Catalog #: D0632



REAGENT

CoA

by [Sigma Aldrich](#)

Catalog #: C4282



REAGENT

Folinic Acid

by [Sigma Aldrich](#)

Catalog #: F7878



REAGENT

Spermidine

by [Sigma Aldrich](#)

Catalog #: 85558



REAGENT

NAD (β -Nicotinamide
adenine dinucleotide
hydrate)

by [Sigma Aldrich](#)

Catalog #: N6522



REAGENT

HEPES

by [Sigma Aldrich](#)

Catalog #: H6147



REAGENT

ATP (Adenosine 5'-
triphosphate dipotassium
salt hydrate)

by [Sigma Aldrich](#)

Catalog #: A8937



REAGENT

GTP (Guanosine 5'-
triphosphate sodium salt
hydrate)

by [Sigma Aldrich](#)

Catalog #: 10106399001 ROCHE

REAGENT
 UTP (Uridine 5'-triphosphate trisodium salt dihydrate)
 by [Sigma Aldrich](#)
 Catalog #: 94370

REAGENT
 Maltodextrin (4.0-7.0 dextrose equivalent)
 by [Sigma Aldrich](#)
 Catalog #: 419672

5 In a 15 mL Falcon tube add the following volumes of the stock solutions prepared before:


Stock solution	Volume needed
HEPES pH 8	1000 µl
Nucleotide mix	396 µl
tRNA solution	160 µl
NAD	76.6 µl
CO-A	160 µl
cAMP	46 µl
Folinic acid	80 µl
spermidine	34 µl
Maltodextrin	2000 µl
TOTAL	3952.6 µl

6 Dissolve well the maltodextrin-based energy solution using vortex. Aliquot in 2 ml Eppendorf tubes and store at -80°C prior to use.

Prepare hexametaphosphate solution

7 Weight 0.15 g of sodium hexametaphosphate and put it into a 15 ml Falcon tube. Dissolve well and fill with water to 5 ml. Put the solution in a boiling water bath for 5 minutes. Kepp in the solution in the lab bench until it reaches room temperature. Aliquot in 1.5 ml Eppendorf tubes and keep at -80°C before use.

REAGENT
 sodium hexametaphosphate
 by [Sigma Aldrich](#)
 Catalog #: 305553

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited