

Carbon isotope discrimination - part I - trapline - v0.1

Stephen Driever, Steven Burgess

Abstract

Protocol for collecting gas for measurement of carbon isotope discrimination during photosynthesis.

See the following references for more background:

H. Griffiths (1993) *Photosynthesis and Production in a Changing Environment* pp 181-1927

Farquhar et al. (1989) *Annu. Rev. Plant Physiol. Plant Mol. Bioi.* 1989. 40:503-37

Citation: Stephen Driever, Steven Burgess Carbon isotope discrimination - part I - trapline - v0.1. **protocols.io**
dx.doi.org/10.17504/protocols.io.qywdxxe

Published: 20 Jun 2018

Protocol

Setting up IRGA

Step 1.

Set the following environmental variables in the LICOR LI-6800

CO2S	400 ppm
Txchg	25 °C
H2OR	10-14 mbar (aim for Tleaf 25 °C)
PAR	2000 (100% Red)
Flow	250 ml min ⁻¹
Fan speed	14500 rpm
deltaP	0.1 kPa

📌 NOTES

Steven Burgess 15 Jun 2018

Control temperature on block (Txchg), not leaf otherwise photosynthetic rate will start to oscillate due to constant adjustment of stomata.

Tleaf should ideally be within +/- 0.03°C of 25°C as measurement is temperature sensitive. 25°C is the value that C3 photosynthetic calculations are parameterized for. You can manipulate with humidity or Txchg on reference

Preparing plant for measurement


Step 2.

- Clamp leaf in cuvette, under a bright ambient light.
- Adjust H2OR so that Tleaf is around 25 °C. This will need adjusting when leaf is reaching stability in the cuvette.

📌 NOTES

Steven Burgess 14 Jun 2018

For tobacco a 4 week old plant prior to flowering is in the ideal state (health example below).



Make sure not to stress the plant by introducing tension when clamping.

Illuminate from above (at least $200\mu\text{E}$), this makes sure the surrounding leaf tissue is photosynthesising to reduce chance of diffusion of metabolites influencing measurements, this is more important for homobaric leaves such as tobacco.

It helps to let the plant acclimate in the lab for $\sim 24\text{h}$ prior to measurement.



Wait with measurements and CO_2 trapping until complete steady state has been reached.

Assimilation rate (A) and steady state fluorescence do not change over time (10 min window)

Stomatal conductance (g_s) does not change over time (10 min window)

This can take $\sim 2\text{h}$ depending on the plant. Check regularly (\sim every 15 mins after the first half hour)

Note that mesophyll conductance can be altered by a whole range of factors (cell size, leaf expansion, age of leaf etc) so care should be taken in performing comparisons and should try to keep to plants at a consistent developmental stage and growth conditions.

Step 3.

When plant has reached a steady state switch inlet gas to 2% oxygen. Use a flow controller to moderate gas use, setting flow to 6.

Change the constant value in LICOR

📌 NOTES

Steven Burgess 15 Jun 2018

2% O₂ is used to suppress photorespiration simplifying downstream calculations of carbon flux. Note this will change assimilation and photosynthesis values which will take time to stabilize again, you may witness oscillations in A and F which will dampen down over time.



When there is no flow the value should be ~0.5 or less



If the cylinder is running low order a replacement, always keep a backup cylinder to avoid running out of gas



Setting up preline

Step 4.

Turn on cold finger

📌 NOTES

Steven Burgess 14 Jun 2018

This will be used to cool down acetone to below -70°C as required for sublimation of frozen CO_2



There is only an ON/OFF setting.

Step 5.

Place cold finger in dewar of acetone and wait to cool down ($>-70^{\circ}\text{C}$), it should turn milky white

📌 NOTES

Steven Burgess 14 Jun 2018




Cover the top of the dewar with foil to prevent ice forming and falling into the acetone. Do not rest the finger on the bottom of the dewar otherwise an ice lump forms and it will stick

Step 6.

Check temperature of acetone is below -70°C

NOTES

Steven Burgess 14 Jun 2018



Step 7.

Fill two dewars 3/4 full with LN₂

NOTES

Steven Burgess 14 Jun 2018



These will be used to condense gases for sample collection. This can be done whilst waiting for the acetone to cool down

Step 8.

Open glass valves and wait until vacuum gauge reads around 5.5 mbar on each prep line. Adjust by opening/closing pump valves if required.

NOTES

Steven Burgess 15 Jun 2018



Provide support from below when opening/closing valve. Apply downward pressure while twisting to minimize the chance of introducing leaks (this is the most common source in the system).



Adjust the vacuum pressure by opening the valve to the pump (red arrow). This ensures that approximately the same amount of gas is collected for reference and sample and ensures consistency between measurements



Step 9.

Check that both collection tubes are open

Step 10.

Check if leaf is still in stable state (both A and gs)

Start Measurement

Step 11.

Match LI-6800 manually

Check settings


- automatch OFF
- Fluorescence OFF

Step 12.

Immerse both traps in LN2 (2 dewars)

+ NOTES

Steven Burgess 15 Jun 2018



The distillation tube should ideally be covered close to the top



Step 13.

Set timer and start logging every minute min for 5 min (LI-6800, autoprogram)

+ NOTES

Steven Burgess 15 Jun 2018



Step 14.

Take fluorescence measurement (F_m' / F_s)


Step 15.

Close the glass valves, watch the vacuum drop to low 10^{-2} or high 10^{-3} E.

NOTES

Steven Burgess 15 Jun 2018

Twist perpendicular to the flow direction otherwise there can be leaks, check the pressure gauge




Step 16.

Close the pump valves.

NOTES

Steven Burgess 15 Jun 2018




Step 17.

For each prep line (one after the other) make sure the collection tubes are open.

NOTES

Steven Burgess 15 Jun 2018



Step 18.

In a single motion remove distillation tube from LN2 and submerge in cold acetone to sublimate gases

NOTES

Steven Burgess 15 Jun 2018



The pressure in the prep line should rise as gas is released

Step 19.

Quickly submerge the collection tube in LN2 to condense the released gas

NOTES

Steven Burgess 15 Jun 2018



The pressure should fall again as gas is drawn down into the collection tube. By following this procedure it is possible to collect gases and remove the water vapour from samples which would interfere with measurements



Collection tube should be dipped in LN2, covering at least the bottom 2cm

Step 20.

Close collection tube


Step 21.

Remove acetone and liquid N dewars

📌 NOTES

Steven Burgess 15 Jun 2018

When the vacuum stabilizes all gas has been condensed again, this is the point when you can close tube and remove from LN2



Collect drips of acetone from distillation tube in beaker



Step 22.

Remove collection tube from prep line and leave to thaw in rack

Put new (empty and open) collection tube in place and open pump valve

🔍 NOTES

Steven Burgess 15 Jun 2018





Opening the pump valve cleans out the prep line and residual air in the new collection tube

Step 23.

Close low O₂ cylinder and adjust constant value in LICOR
