



Jun 09,  
2019

## Cleaning up a biohazardous spill inside a centrifuge

In 1 collection

Steven Wilhelm<sup>1</sup>, Gary LeClerc<sup>1</sup>, Ashley Humphrey<sup>2</sup>

<sup>1</sup>The University of Tennessee, Knoxville, <sup>2</sup>University of Tennessee, Knoxville

Working dx.doi.org/10.17504/protocols.io.wgkfbuw

The Aquatic Microbial Ecology Research Group - AMERG (The Buchan, Zinser and Wilhelm labs)



Ashley A. Humphrey  
University of Tennessee, Knoxville



### ABSTRACT

**BIOLOGICAL SPILL RESPONSE** method to be used in the event of a biosafety incident.



SpillCleanupProcedure\_20  
17.07.27 (1).pdf

### SAFETY WARNINGS

For chemical safety, please refer to the SDS (Safety Data Sheet).

- 1 Close off the area and allow aerosols to settle.
- 2 Notify others including supervisor.
- 3 Wait ⌚ 00:30:00 to allow the aerosols to settle.
- 4 Don appropriate PPE: laboratory coat, safety glasses and Nitrile gloves.
- 5 Remove rotors and bucket and place in Biosafety Cabinet.
- 6 Thoroughly disinfect the inside and outside of the centrifuge rotor, cups and accessories and allow proper contact time. After disinfection, move to sink for a thorough rinse, dry thoroughly.
- 7 For the centrifuge interior, follow the steps above for cleaning up a spill outside a Biosafety Cabinet.

- 8 Segregate contaminated cleanup materials into the appropriate biohazardous waste containers.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited