

Restriction digest

Snehadri Sinha

Abstract

Citation: Snehadri Sinha Restriction digest. **protocols.io**

dx.doi.org/10.17504/protocols.io.gzmbx46

Published: 17 Jan 2017

Protocol

Step 1.

Combine the materials:

>500 ng DNA (500 ng for diagnostic digest, >1000ng for restriction cloning)

0.5-1 uL each restriction enzyme

Buffer (appropriate buffer indicated by enzyme manufacturer, to a final concentration of 1x)

dH₂O to total volume of 20 ul.

Step 2.

Mix gently by pipetting.

Step 3.

Incubate tube for an appropriate temperature at an appropriate time (usually 37 C, time varies; generally 1 h for NEB enzymes, 10-30 min for FastDigest enzymes. Follow the manufacturer's instructions.).

Step 4.

Visualize the results of your digest, conduct gel electrophoresis.