

Capped RNA Synthesis (E2050)

New England Biolabs

Abstract

The kit formulation allows for efficient capped RNA synthesis using cap analog (ARCA).

Citation: New England Biolabs Capped RNA Synthesis (E2050). **protocols.io**

dx.doi.org/10.17504/protocols.io.cryv7v

Published: 04 Feb 2015

Guidelines

We strongly recommend wearing gloves and using nuclease-free tubes and reagents to avoid RNase contamination. Reactions are typically 20 µl but can be scaled up as needed. Reactions should be assembled in nuclease-free microfuge tubes or PCR strip tubes.

The kit formulation allows for efficient capped RNA synthesis using cap analog (ARCA). The recommended ratio of cap analog to GTP is 4:1. Increasing the ratio of cap analog to GTP will increase the proportion of capped RNA transcripts; however, it also significantly decreases the yield of the reaction. Cap analogs are sold separately. Please refer to the [companion products](#) section.

Cap analog (ARCA, [NEB #S1411](#)) is supplied in a lyophilized form of 1 µmol per tube. Dissolving it in 25 µl nuclease-free water will yield a concentration of 40 mM.

The yield per reaction is 30–40 µg RNA with approximately 80% capped RNA transcripts. Figure 1 shows the time course of capped RNA synthesis from 1 µg control template. Most reactions will be complete in 1 hour.





Figure 1. Capped RNA Synthesis with ARCA

Reactions were incubated at 37°C in a thermocycler. Transcripts were purified by spin columns and quantified on a NanoDrop Spectrophotometer.

Before start

We strongly recommend wearing gloves and using nuclease-free tubes and reagents to avoid RNase contamination. Reactions are typically 20 µl but can be scaled up as needed. Reactions should be assembled in nuclease-free microfuge tubes or PCR strip tubes.

Materials

 HiScribe T7 Quick High Yield RNA Synthesis Kit - 50 rxns [E2050S](#) by [New England Biolabs](#)

Protocol

Step 1.

Prepare 40 mM cap analog



REAGENTS

 3'-O-Me-m7G(5')ppp(5')G RNA Cap Structure Analog - 1 µmol [S1411S](#) by [New England Biolabs](#)



NOTES

New England Biolabs 04 Feb 2015

Cap analog (ARCA, [NEB #S1411](#)) is supplied in a lyophilized form of 1 µmol per tube. Dissolving it in 25 µl nuclease-free water will yield a concentration of 40 mM.

Step 2.

Thaw the necessary kit components, mix and pulse-spin in a microfuge to collect solutions to the bottoms of tubes

Step 3.

Assemble the reaction at room temperature in the following order:



PROTOCOL

. [Capped RNA E2050 Mixture](#)

CONTACT: [New England Biolabs](#)

Step 3.1.

Nuclease-free water so that total volume is **20 µl**

Step 3.2.

NTP Buffer Mix, **2 µl** (2 mM each NTP final)



AMOUNT

2 µl Additional info:

Step 3.3.

Cap Analog (40 mM), **4 µl** (8 mM final)



AMOUNT

4 µl Additional info:

Step 3.4.

Template DNA, X µl (1 µg)

Step 3.5.

T7 RNA Polymerase Mix, **2 µl**

 **AMOUNT**

2 µl Additional info:

Step 4.

Mix thoroughly and pulse-spin

Step 5.

Incubate at 37°C for 2 hours

 **DURATION**

02:00:00

 **NOTES**

New England Biolabs 04 Feb 2015

The yield per reaction is 30–40 µg RNA with approximately 80% capped RNA transcripts. See guidelines for figure showing the time course of capped RNA synthesis from 1 µg control template.

New England Biolabs 04 Feb 2015

Most reactions will be complete in 1 hour.

Step 6.

Optional step: To remove template DNA, add **2 µl** of DNase I (RNase-free), mix and incubate at 37°C for 15 minutes.

Step 7.

Proceed with purification of synthesized RNA or analysis of transcription products by gel electrophoresis