

How to culture ichthyosporean *Abeoforma whisleri*, a unicellular relative to animals

Maria Rubio-Brotons

Abstract

Citation: Maria Rubio-Brotons How to culture ichthyosporean *Abeoforma whisleri*, a unicellular relative to animals. **protocols.io**

dx.doi.org/10.17504/protocols.io.gqkbvuw

Published: 12 Dec 2016

Protocol

Add 8mL of marine broth into sterile flask

Step 1.

To perform under hood

Marine Broth medium (Ref. Marine Broth 2216, BD-Difco)

Tissue Culture Flask 50 mL capacity (Ref. 353108 Blue Vented Cap, Falcon)

Add 100uL of cells into flask

Step 2.

Passage of 1 week old cells

Mix culture carefully (do not to touch cap)

Step 3.

📌 NOTES

Maria Rubio-Brotons 12 Dec 2016

For optimal culture conditions, perform passages weekly

Store in the incubator at 23°C

Step 4.