



RT-PCR for NoV

PLOS One

Hengyun Guan¹, Chunrong Wang¹, Lanzheng Liu¹, Guoliang Yang¹

¹Jinan Center for Disease Control and Prevention, Jinan, China

dx.doi.org/10.17504/protocols.io.q3sdyne

Hengyun Guan

ABSTRACT

This assay amplified the partial VP1 gene (Region C) of the norovirus genome. The yield products could be sequenced.

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0209245>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Liu L, Guan H, Zhang Y, Wang C, Yang G, Ruan S, Zhao H, Han X (2018) The prevalence of non-GII.4 norovirus genotypes in acute gastroenteritis outbreaks in Jinan, China. PLoS ONE 13(12): e0209245. doi: [10.1371/journal.pone.0209245](https://doi.org/10.1371/journal.pone.0209245)

PROTOCOL STATUS

Working

Reagents

- 1 SuperScript™ III One-Step RT-PCR with Platinum™ Taq
By [Life Technologies](#)
Catalog #: [12574026](#)

RNA extraction

- 2 Extract total RNA of specimens with the MagNA Pure LC Total Nucleic Acid Isolation Kit (Roche, Mannheim, Germany) according to the manufacturer's instructions.

Reaction system

- 3 The primers (G1-SKF/G1-SKR for GI, COG2F/G2SKR for GII) were derived from the studies of Kojima, S. and Shinohara, M. et al.

The reagent composition	Volume25 ul
2×Reaction Mix	12.5 ul
Forward Primer (25 uM)	0.5 ul
Reverse Primer (25 uM)	0.5ul
Taq mix	1.0 ul
Rnase Inhibitor (30-40U/ μl)	0.5 ul
RNA	3 ul
DEPC H ₂ O	to 25ul

Amplification

- 4 The amplification conditions were as follows: 42° for 30 min; 95° for 15 min, 40 cycles of 95° for 30 s, 50° for 30 s, and 72° for 30 s and a final extension of 72° for 10 min.

Products

- 5 The expected amplicon size of GI was 330 bp and GII was 387 bp. The yield products were sequenced to validate the the genotype of NoV.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited