



May 15,
2019

Working

UC Davis - Metabolomics: Sample preparation for GCTOF analysis [↗](#)

Oliver Fiehn¹

¹University of California, Davis

[dx.doi.org/10.17504/protocols.io.ytjfwkn](https://doi.org/10.17504/protocols.io.ytjfwkn)

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org

Lili Liang

ABSTRACT

Summary:

This SOP describes sample extraction and sample preparation for primary metabolism profiling by gas chromatography / time of flight mass spectrometry (GCTOF)

References:




Fiehn O, Kind T (2006) Metabolite profiling in blood plasma. In: Metabolomics: Methods and Protocols. Weckwerth W (ed.), Humana Press, Totowa NJ

EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=121&docType=Protocol>

MATERIALS

NAME	CATALOG #	VENDOR
Centrifuge	5415D	Eppendorf Centrifuge
Calibrated pipettes 1-200 ul and 100-1000 ul	View	
Eppendorf tubes 1.5 mL uncolored	022363204	Eppendorf Centrifuge
ThermoElectron Neslab RTE 740 cooling bath	RTE 740 cooling bath	
MiniVortexer	58816-121	VWR Scientific
Orbital Mixing Chilling/Heating Plate		Torrey Pines Scientific Instruments
Speed vacuum concentration system		Labconco Centrивap cold trap
Precision balance with accuracy ± 0.1 mg		
2mL crimp vials with Target Micro-Serts		
Agilent Electronic crimper and decapper	View	Agilent Technologies
Acetonitrile LCMS quality	9829-02	JT Baker
Isopropanol HPLC solvent	9095-02	JT Baker
pH paper 5-10	108027	EMD Chem. Inc.
Nitrogen line with pipette tip		
Methoxyamine hydrochloride [MeOX]	226904	Sigma Aldrich
Pyridine	270970-4X25ML	Acros Organics

NAME 	CATALOG # 	VENDOR 
N-methyl-N-(trimethylsilyl)-trifluoroacetamide [MSTFA]	394866	Sigma Aldrich
FAME markers (refer to FAME marker SOP for preparation)		

MATERIALS TEXT

Note:

Eppendorf, [RRID:SCR_000786](#)

Sigma-Aldrich, [RRID:SCR_008988](#)

BEFORE STARTING

Starting material:

Plasma/serum: 30 µl sample volume or aliquot

1 Preparation of extraction mix before experiment:

- (1). Check pH of acetonitrile and isopropanol (pH7) using wetted pH paper
- (2). Acetonitrile, isopropanol and water are mixed in volumes in proportion 3 : 3 : 2
- (3). Rinse the extraction solution mix for 5 min with nitrogen with small bubbles. Make sure that the nitrogen line was flushed out of air before using it for degassing the extraction solvent solution

2 Sample Preparation:

- (1). Switch on bath to pre-cool at -20°C ($\pm 2^{\circ}\text{C}$ validity temperature range)
- (2). Gently rotate or aspirate the blood samples for about 10s to obtain a homogenised sample.
- (3). Aliquot 30µl of plasma sample to a 1.0 mL extraction solution. The extraction solution has to be pre-chilled using the ThermoElectron Neslab RTE 740 cooling bath set to -20°C .
- (4). Vortex the sample for about 10s and shake for 5 min at 4°C using the Orbital Mixing Chilling/Heating Plate. If you are using more than one sample, keep the rest of the sample on ice (chilled at $<0^{\circ}\text{C}$ with sodium chloride).
- (5). Centrifuge samples for 2min at 14000 rcf using the centrifuge Eppendorf 5415 D.
- (6). Aliquot two 450µL portions of the supernatant. One for analysis and one for a backup sample. Store the backup aliquot in -20°C freezer.
- (7). Evaporate one 450µL aliquots of the sample in the Labconco Centrивap cold trap concentrator to complete dryness.
- (8). The dried aliquot is then re-suspended with 450 µL 50% acetonitrile (degassed as given above).
- (9). Centrifuged for 2 min at 14000 rcf using the centrifuge Eppendorf 5415.
- (10). Remove supernatant to a new Eppendorf tube.
- (11). Evaporate the supernatant to dryness in the Labconco Centrивap cold trap concentrator.
- (12). Submit to derivatization.

3 Derivatization

- Prepare 40mg/mL MeOX solution in pyridine. Weigh out methoxyamine hydrochloride in 1.5mL Eppendorf tube on balance and add appropriate amount of pyridine.

- Vortex MeOX solution and sonicate at 60°C for 15 minutes to dissolve.
- Add 10 µL of 40mg/mL Methoxyamine hydrochloride [MeOX] solution to each dried sample and standard
- Shake at maximum speed at 30°C for 1.5 hours.
- To 1mL of MSTFA and add 10 µL of FAME marker. Vortex for 10sec.
- Add 91 µL of MSTFA + FAME mixture to each sample and standard. Cap immediately.
- Shake at maximum speed at 37 °C
- Transfer contents to glass vials with micro-serts inserted and cap immediately.
- Submit to GCTOF MS analysis

4 Quality assurance

- For each sequence of sample extractions, perform one blank negative control extraction by applying the total procedure (i.e. all materials and plastic ware) without biological sample.
- Use one commercial plasma/serum pool sample per 10 authentic subject samples as control. If no suitable commercial blood sample is available, prepare a large pool sample during the thawing/mixing step by aliquoting 100 ul per 1 ml plasma sample, and aliquot such pool sample for 1 pool extract per 10 authentic subject samples.
- Prepare at least one NIST plasma extract in the same manner

***DISPOSAL OF WASTE: Collect all chemicals in appropriate bottles and follow the disposal rules.
Collect residual plasma / serum samples in specifically designed red 'biohazard' waste bags.***



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited