

# Simple electroporation for efficient CRISPR/Cas9 genome editing in murine zygotes

Simon E. Tröder, Branko Zevnik

## Abstract

Electroporation of zygotes represents a rapid alternative to the elaborate pronuclear injection procedure for CRISPR/Cas9-mediated genome editing in mice. However, current protocols for electroporation either require the investment in specialized electroporators or corrosive pre-treatment of zygotes which compromises embryo viability. Here, we describe an easily adaptable approach for the introduction of specific mutations in C57BL/6N mice by electroporation of intact zygotes using a common electroporator with synthetic CRISPR/Cas9 components and minimal technical requirement.

**Citation:** Simon E. Tröder, Branko Zevnik Simple electroporation for efficient CRISPR/Cas9 genome editing in murine zygotes. **protocols.io**

dx.doi.org/10.17504/protocols.io.ndzda76

**Published:** 08 Mar 2018

## Guidelines

Use only embryo-grade reagents.

M2 and M16 media are prepared as published (Behringer, R., et al., (2014) Manipulating the mouse embryo : a laboratory manual, Fourth edition. ed. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York).

Commercial M2 (Sigma, #M7167) and KSOM/AA (Merck Millipore, #MR-106-D) can be used instead.

Mouse pre-implantation embryos are incubated in at least 4h pre-equilibrated M16 or KSOM/AA in a CO<sub>2</sub> incubator (5% CO<sub>2</sub>, 37°C, 95% humidity).

Embryo-grade T<sub>10</sub>E<sub>0.1</sub> (10 mM Tris-HCl, 0.1 mM EDTA, pH 7.4) buffer is prepared as described (Chu, V.T., et al., (2016) Efficient generation of Rosa26 knock-in mice using CRISPR/Cas9 in C57BL/6 zygotes. BMC Biotechnol 16, 4.).

Mouse zygotes are collected from the oviducts of superovulated females as described in published protocols (Behringer, R., et al., (2014) Manipulating the mouse embryo : a laboratory manual, Fourth edition. ed. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York).

## Materials

- Cas9 Nuclease 1074181 by [Integrated DNA Technologies](#)
- tracrRNA 1072532 by [Integrated DNA Technologies](#)
- crRNA by [Integrated DNA Technologies](#)
- ssODN (Ultramer DNA Oligonucleotides) by [Integrated DNA Technologies](#)
- Opti-MEM (Reduced Serum Medium) 31985062 by [Thermo Fisher Scientific](#)
- ✓ M2 Medium by Contributed by users
- ✓ M16 Medium by Contributed by users
- ✓ T10E0.1 Buffer by Contributed by users
- Electroporation Cuvette 1mm 1652089 by [BioRad Sciences](#)
- 60 mm Center Well Organ Culture Dish 353037 by [Corning](#)

## Protocol

### guide RNA annealing

#### Step 1.

Resuspend lyophilized crRNA, tracrRNA and ssODN in T<sub>10</sub>E<sub>0.1</sub> buffer to 100 µM (e.g. 5 nmol in 50 µl)

(Store at -80 °C until use)

### guide RNA annealing

#### Step 2.

Combine 5 µl crRNA (100 µM) and 5 µl tracrRNA (100 µM) in a nuclease-free PCR tube to yield an equimolar crRNA:tracrRNA duplex solution of 50 µM

### guide RNA annealing

#### Step 3.

Heat to 95 °C for 5 min and cool down at 5 °C/ min in a thermocycler

(crRNA:tracrRNA duplex can be stored for weeks at -80 °C)

### Preparation of the electroporation mix

#### Step 4.

Add 1.6 µl crRNA:tracrRNA duplex and 1.3 µl Cas9 nuclease to 15.1 µl Opti-MEM in a nuclease-free tube and vortex

### Preparation of the electroporation mix

#### Step 5.

Incubate mix at room temperature for 10 min

#### Preparation of the electroporation mix

##### Step 6.

Place tube on ice, add 2 µl ssODN and vortex

#### Preparation of the electroporation mix

##### Step 7.

Quick-spin at 4 °C and keep tube on ice until use

#### Preparation of the electroporation mix

##### Step 8.

Summary of the 20 µl electroporation mix:

Reagent	Stock concentration	Final concentration	Volume
crRNA:tracrRNA duplex	50 µM	4 µM	1.6 µl
Cas9 nuclease	61 µM (10 µg/µl)	4 µM	1.3 µl
ssODN	100 µM	10 µM	2.0 µl
Opti-MEM	-	-	15.1 µl

#### Electroporation of zygotes

##### Step 9.

Collect zygotes from the oviducts of superovulated females as described in published protocols

#### Electroporation of zygotes

##### Step 10.

Wash the zygotes in five drops of M2

#### Electroporation of zygotes

##### Step 11.

Wash up to 50 zygotes in one drop of Opti-MEM

#### Electroporation of zygotes

##### Step 12.

Transfer zygotes with as little media as possible to the 20 µl electroporation mix

#### Electroporation of zygotes

##### Step 13.

Using a 20 µl pipette transfer the entire drop including the zygotes into a pre-warmed (37 °C) 1 mm electroporation cuvette

(Ensure retrieving all zygotes by quickly aspirating the entire drop. Slow aspiration will leave zygotes behind)

#### Electroporation of zygotes

##### **Step 14.**

Insert cuvette into a standard electroporator (e.g. BioRad Gene Pulser Xcell electroporator)

#### Electroporation of zygotes

##### **Step 15.**

Apply two square wave pulses at 30 V and 3 ms duration with a 100 ms interval

#### Electroporation of zygotes

##### **Step 16.**

Retrieve the zygotes by flushing the cuvette with 100 µl pre-incubated M16 using a 100 µl pipette into a culture dish (e.g. 60 mm Center Well Organ Culture Dish)

#### Electroporation of zygotes

##### **Step 17.**

Wash the cuvette with 100 µl pre-incubated M16

#### Electroporation of zygotes

##### **Step 18.**

Transfer all zygotes to a new culture dish with 500 µl pre-incubated M16

#### Electroporation of zygotes

##### **Step 19.**

Incubate zygotes in M16 until the two cell stage and transfer the developed embryos into pseudopregnant foster mice

(Embryos may also be transferred at the one cell stage)